

Product datasheet for **MC200390**

Xlr (NM_011725) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Xlr (NM_011725) Mouse Untagged Clone
Tag: Tag Free
Symbol: Xlr
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >BC005560 sequence for NM_011725

```
GTGGCCACGAAGCCGAGTTCTTGAGGAAGCGAAGTTCTTGAGGAAGCCGAGAAGTTCTTGAGAAAACCTA  
AATTTCTTAGTGGAAAGGCAGTGTCTTGAATATTCTCTAGAGACAACAATGGAAAACCTGGGACTTGTC  
GTGATGAAATGCAAGATGGGAATGCTCCAGAATTGGACGTTATTGAAGAACATAATCCAGTAACTCGTGA  
TGATGAGAATGCAATCCTGAAGAAGTAGTTGGAGATACACGATCTCCAGTACAAAATATTCTGGGAAAA  
TTTGAAGGTGACATTAACAAGAGGCTTCACATAAAGAGAAAACGCATGGAAAACCTATATCAAGATTCTT  
TCAAAGACAGCAACGTGAAATTAGAACAACCTTTGGAAAACGAACAACAAGAGAGAAAGAAGATTAACAA  
CAAATTTGTGAGCAGTATATACTACATTTTCAAGAGTTGATATGGATGTCCAGAAATTCATGAAGAA  
CAAGAAAAATCAGTGAACAATTATCAAAAAGAACAACAAGCATTAAAACGTCCAAATGTAGTCAGAGCC  
AGACCTGGAAGCAATTAAGACATGCATGAAAACCTACATGGAGGGTCTGATGAACTTGGAGACCAACAA  
CTACAATATGCTTTTGTAGATGGTGAACCTGAGAAAAGAAATGTCAGTGTAAAAAAGACCTCATG  
AAGCATACTCTGAAGTACTCTTCTAGTTTCCCATCTTCAGACTAAGGAAGCATGATATTTGCACCTTGC  
TGGTACATTTGTAACCAATAAAAAATTCGCTAACTTAAAAA
```

Restriction Sites: RsrII-NotI
ACCN: NM_011725
Insert Size: 627 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC005560](#), [AAH05560](#)

RefSeq Size: 822 bp

RefSeq ORF: 627 bp

Locus ID: 22441

UniProt ID: [P05531](#)

Cytogenetics: X A5