

Product datasheet for **MC200431**

Tspan4 (NM_053082) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tspan4 (NM_053082) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Tspan4
Synonyms:	AI325509; AI746565; D130042I01Rik; NAG-2; Tm4sf7; Tspan-4
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC003482 sequence for NM_053082 ATCAGGCTCAAGATCTGGACTGAAGACTGCTGGTGGCCAGAGAGCTGAAGCGCTGCAGCATGGCGCGCGG CTGCCTCCAGGGCGTCAAGTACCTCATGTTCCGCTTCAATCTGCTCTTCTGGCTGGGTGGCTGTGGTGT CTGGGTGTTGGCATCTGTTGGCTGCCACACAGGAACTTTGCCACCTTATCATCCTCATTCCATCCT TGTCGGCTGCCAACCTGCTCATCGTCACCGGGACCTTCGTCATGGCCATCGGCTTCGTGGCTGCATTGG GGCCCTCAAGGAGAACAAGTGCCTACTGCTCACTTTCTTTGTGCTGCTGCTAGTGTTCCTGCTGGAA GCCACCATTGCTGTGCTCTTCTTTGCCTACAGTGACAAGATTGACAGTTATGCCCAACAAGACCTGAAGA AGGGCCTGCATCTGTATGGCACACAGGGCAACGTGGGCCACCAATGCCTGGAGCATCATCCAGACTGA TTTCCGATGCTGTGGAGTTTCCAATTACACTGATTGGTTTGAGGTATAACAATGCCACTCGTGTGCCTGAC TCCTGCTGTCTGGAGTTCAGTGATAGCTGTGGTTACATGAACCTGGTACCTGGTGAAGTCGCCCTGTT ATGAGACAGTGAAGGCCCTGGCTCCAGGAGAACCTGCTAGCTGTGGGCATCTTTGGACTGTGCACGGCACT GGTGCAGATTCTGGGCCCTCACCTTCGCTATGACCATGTAAGTGCAGGTAAAGGCGGACACCTACTGT GCATAGGTCTCACAGCCTCTGCTTCTCTGCCAAAGGACCCGCCACGGGAAGGAAGATGTCCAAACCCGC TCCCCCTGCAGGACTGGTGTCTGCTGAATGCAGGGGGAAGAGCTGGAGGCCAGAGCACCAGCATCCAT CCCTGCACCTGAGGGCACTCTGCTTGGGTGGCAGCAGACACCAGACTATAGAGTTGGTGCCAAATGGCAC CCAGGGCTAAGGGCCAGGAACTTTGACACTTTTATATCCATATATACCTTCCAGTGTTCACACAGGA TGGGTTCCAGCCTTGCATGGTTGCTAGGTTGGGGCAGGCCAGTCTGCTCAAGGCCACGATGATCT TGGTGTCCAGGCAGGGCCAGGTCCCATGATCCACATTCCATGACAGTGTCTGGTTCAGATGCATCTT AATAAAGTGTGTGAGCAGACTGTGGTCTTGCCAACACACAAAGTTAAAAAAAAAAAAAAAA
Restriction Sites:	RsrII-NotI
ACCN:	NM_053082
Insert Size:	717 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC003482 , AAH03482
RefSeq Size:	1252 bp
RefSeq ORF:	717 bp
Locus ID:	64540
UniProt ID:	Q9DCK3
Cytogenetics:	7 F5