

Product datasheet for **MC200452**

Tspan1 (NM_133681) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tspan1 (NM_133681) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Tspan1
Synonyms:	9030418M05Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC003448 sequence for NM_133681 CCTGGAGCTGTTCTGTTCCCTTCTGGAGTGCAACAGATCCCCTTCCAGAACCTATTGCCAAGAGCCCTGAAG CAGGACACACCATGCAGTGTTCAAATTCATTAAGGTCATGATGTTCTTCAATCTACTCATCTTTCT CTGTGGTGCAGCCCTGTTGGCTGTGGGAATCTGGGTGTCCGTCGATGGGACATCTTCTGAAGGCTTC GGATCACTATCATCCAGTCCATGCAGTTGTCAACGTGGGCTACTTCTCATCGCCGCTGGTGCTGTGC TCTTATTCTTGGTTTCTGGGCTGCTATGGTGTCTACTCTGAGAACAAGTGTGTGCTCATGATGTTCTT TTCCATCCTCCTCATCATCTTTCATCGTGAGATTGCAGGTGCTGTGGTTGCTTTGGTGTACACCACATTG GCTGAACAATTCCTGACTCCTGGTGGTGCCTGCTATCGAAAAAGACTATGGTTACCAGACTGATTTCA CCCAAGTATGGAACACTACAATGGAAGAGTTGCATTGCTGTGGCTTTAACAACTACACAGATTTTAAATGC CTCACGTTTCGTCAAAGAGAATAAAGTCTTCCCCCACCCTGTTGTGCCAACCTGGCAACCATACAGTT GAACCATGCACCGAGGAGAAGGCCAAAAGTATGAAAGTACAGGGTGTGTTCAAAGAGATTCTGCATAGAA TCAGAGCCAATGCAGTCACTGTGGGTGGTGTGGCAGTTGGAGTTGCGGCCCTAGAGCTGGCTGCCATGGT TGTATCCATGTATCTATACTGCAATCTGAAATAAGACTACTACTTCTCCTGACTTGCTGCCGTCACATG GAAACCAGGAAGAGGCACTTGGCACAGCTAAGCAGCAGCAATGGGTGGAGATGAAATGTCATTTGGGCTG GAATGGATTTGCCTTTGCTGCTCAGGGACTTTGGGCTAAATAGGGACCACCTCGTGGGCTTTGCTGCTCC TGAGACTCTTGATGATCCTAGCACTCCACTAGTGGATGAAGGCACTTACAGCCTGGGCATTTTCTAACAC TGTAAGTCAATGAAGGAATGTATAGAAGGCATTTTGGAACTATAAACCAATGAAAATATTTCTAATTG ACCTGTAAAAAAAAAAAAA
Restriction Sites:	RsrII-NotI
ACCN:	NM_133681
Insert Size:	723 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC003448</u> , <u>AAH03448</u>
RefSeq Size:	1209 bp
RefSeq ORF:	723 bp
Locus ID:	66805
UniProt ID:	<u>Q99J59</u>
Cytogenetics:	4 D1