

Product datasheet for **MC200501**

Dyrk3 (NM_145508) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dyrk3 (NM_145508) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Dyrk3
Synonyms:	BC006704
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >BC006704 sequence for NM_145508
 AGCGAAGCCAGCAGGTGCGGTGCTGGAGGCGACCCACCGGGCCAACCGCCGCTCTGTGTCCGCGGGGCG
 AGCTCGCGATGGGAGGCGCAGCCCGCATCGCGGGAGGAAGGACGCGGCTGCCGGGGCCGGGCTCCC
 GCCGCAGCAGCGGAGGTTGGGGGATGGTGTCTATGATACTTTTCATGATGATAGATGAAACCAAGGGCCCA
 CCCTATTCGGACACATTCAGCAACCCCTCTGAAGCACCTGTCTCCAGAAGGCTAAATATTACCCTGAGC
 CACTCACGAGAGGCCACACTCAGCACTTTGTGAATGGAAGTGAGATGAAGGTGGAACAACCTTTTCAAGA
 ATTTGGCAACAGAAGATCCAATACTCTCCAGTCAGATGGCATCAGCAACTCGGAAAAGTCTTCTCCGGCT
 TCTCAGGGGAAAAGTTTCAGAAAAGCCTGAGCGCGGTGAAGTGCAACCTTTCTCCAGGCCGTCTAAGGTGC
 TGCCGCTGACTCCTGAGCAAGCGCTGAAGCAGTATAAACACCCACCTCACTGCCTACGAGAAGCTGGAGAT
 CGTCAGCTACCCAGAAATCTACTTTGTGGGTCCGAATGCCAAAAAGCGGCAGGGAGTTATTGGTGGTCCC
 AATAACGGTGGCTACGACGATGCGGATGGGGCCTATATTCATGTGCCCGGGACCATCTGGCTTACCCT
 ATGAGGTGCTGAAAATCATCGGCAAGGGGAGTTTTGGACAGGTAGCCCGGTCTATGATCACAACTTCG
 GCAGTACGTGGCCCTGAAAATGGTGCAGTATGAGAAACGTTCCATCGCCAGGCAGCCGAGGAGATCCGG
 ATCTTGGAGCATCTTAAAAAGCAAGACAAAACCTGGTAGCATGAACGTATCCACATGCTAGAAAATTTCA
 CCTTCCGGAACCGTGTGCATGGCCTTTGAACTGCTAAGCATAGACCTGTATGAGCTTATTAACAAAAA
 CAAGTTCAGGGCTTCAGCGTCCAGTTGGTTCGGAAGTTCGCCAGTCCATCCTACAGTCCCTGGACGCT
 CTGCACAAAAACAAGATCATCCACTGCGACCTGAAGCCGAAAACATTCTCTGAAACACCACGGCCGAA
 GCGCCACCAAGGTCATCGACTTTGGCTCCAGCTGCTTCGAGTATCAGAAGCTTTACACGTATATCCAGTC
 CCGCTTACAGAGCCCCAGAGATCATCTGGGGTGGCGCTACAGCACCCCGATCGACATATGGAGCTTT
 GGTTCATCCTTGCAAACTTTTACAGGACAGCCCTGTTCCCTGGAGAGGACGAAGGAGACCAGTTGG
 CCTGCATGATAGAGTTGCTAGGAATGCCACCGCAGAACTTCTGGAGCAATCCAAGCGTCCAAAGTACTT
 TATTAACCTCAAAGGCTTGCCCTGATACTGCTCCGTATCTACCCAGACGGACGGGAGGGTGGTCTTCTC
 GGGGTCGCTCACGCAGGGGTAAGAAAGCGAGCCCGCCAGGCAGCAAAAGACTGGGCAACCGCACTGAAGG
 GCTGTGGTGACTACTTGTTCATAGAGTTTCTGAAACGATGCCTCCAGTGGGACCCCTCTGCCCGCTCAC
 CCCGGCTCAAGCATTAAAGACATCCTTGGATTAGCAAGTCTACACCCAAACCTCTACCATGGACAAGGTG
 CCAGGGAAGCGGGTAGTTAACCTACAAATGCTTTCAGGGACTGGGTTCCAAGTGCCTCCAGTCGTTG
 GGATAGCCAGTAAGCTTAAAGCTAACCTAATGTCCGAAACAGTGGTAGTATACCTCTGTGCAGTGTATT
 GCCAAAGTGTATTAGCTAGTGGACCACTCAGAGACTGATACATATCATATGATTTTTTAATTACCTTGCA
 AACATGCAAAATGGAAAACGGAATAATTGAAGCCATTCACTGATGGATATGTTTTTGTAGACTTTTTTT
 TAACAAGGCAGAACATTTTTATGACTATAAAAGAACGCTTCAAGGGCTAATGTCAAACAGCTGTAT
 TGCCATCTGGAGTATACATTAATGACTTTTTTCATAGGTCAAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: RsrII-NotI

ACCN: NM_145508

Insert Size: 1761 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC006704](#), [AAH06704](#)

RefSeq Size: 2091 bp

RefSeq ORF: 1761 bp

Locus ID: 226419

UniProt ID: [Q922Y0](#)

Cytogenetics: 1 E4

Gene Summary: Dual-specificity protein kinase that promotes disassembly of several types of membraneless organelles during mitosis, such as stress granules, nuclear speckles and pericentriolar material (By similarity). Dual-specificity tyrosine-regulated kinases (DYRKs) autophosphorylate a critical tyrosine residue in their activation loop and phosphorylate their substrate on serine and threonine residues (PubMed:12356771). Acts as a central dissolvase of membraneless organelles during the G2-to-M transition, after the nuclear-envelope breakdown: acts by mediating phosphorylation of multiple serine and threonine residues in unstructured domains of proteins, such as SRRM1 and PCM1 (By similarity). Does not mediate disassembly of all membraneless organelles: disassembly of P-body and nucleolus is not regulated by DYRK3 (By similarity). Dissolution of membraneless organelles at the onset of mitosis is also required to release mitotic regulators, such as ZNF207, from liquid-unmixed organelles where they are sequestered and keep them dissolved during mitosis (By similarity). Regulates mTORC1 by mediating the dissolution of stress granules: during stressful conditions, DYRK3 partitions from the cytosol to the stress granule, together with mTORC1 components, which prevents mTORC1 signaling (By similarity). When stress signals are gone, the kinase activity of DYRK3 is required for the dissolution of stress granule and mTORC1 relocation to the cytosol: acts by mediating the phosphorylation of the mTORC1 inhibitor AKT1S1, allowing full reactivation of mTORC1 signaling (By similarity). Also acts as a negative regulator of EPO-dependent erythropoiesis: may place an upper limit on red cell production during stress erythropoiesis (By similarity). Inhibits cell death due to cytokine withdrawal in hematopoietic progenitor cells (By similarity). Promotes cell survival upon genotoxic stress through phosphorylation of SIRT1: this in turn inhibits p53/TP53 activity and apoptosis (PubMed:20167603).[UniProtKB/Swiss-Prot Function]