

Product datasheet for **MC200541**

Cpsf3l (BC008240) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cpsf3l (BC008240) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Cpsf3l
Synonyms:	2410006F12Rik; AU044035
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC008240
ACCTGAAGGTGCCATATACTTCTCTACGGCCTGACAGAGAAAGCCAACCACTACTACAAGCTCTTCAT
CACCTGGACCAACCAGAAGATCCGGAAGACGTTTCGTCCAGAGGAACATGTTTGAGTTTAAGCACATCAAA
GCCTTTGACCGGACATTTGCTGACAACCCAGGTCCCATGGTTGTGTTTGCCACACCTGGGATGCTGCATG
CTGGCCAATCCCTGCAAATCTTCCGGAAGTGGCAGGAAATGAGAAGAACATGGTCATCATGCCTGGTTA
TTGTGTGCAAGGCACCGTGGCCACAAGATCCTTAGTGGGCAACGTAAACTGAAATGGAAGGGCGCCAG
ATGCTGGAGGTAAGATGCAAGTGAATACATGTCCTTCAGCGCCATGCTGATGCCAAGGGCATCATGC
AGCTGGTAGGACAAGCGGAGCCGAGAGTGTGCTGTTGGTGCATGGCGAAGCCAAGAAGATGGAGTTCCT
GAGGCAGAAGATTGAGCAGGAATCCGGGTCAGCTGCTACATGCCGGCCAATGGTGAACAGTGACGCTG
CCCACAAGCCCTAGCATCCCTGTGGGCATCTCCCTGGGATTGCTGAAGCGGGAGATGGTGCAGGGGCTGT
TGCCTGAAGCCAAAAAACCTAGGCTTTTACATGGCACCCCTCATTATGAAAGACAGTAATTTCCGGCTGGT
ATCCTCAGAGCAAGCCCTGAAGGAGCTGGGCCTGGCAGAGCACCAACTTCGCTTACCTGTCGTGTGCAC
CTGCAGGACACACGCAAAGAACAAGAAACAGCTTTACGTGTGTACAGCCATCTCAAGAGCACCCCTCAAGG
ACCACTGTGTACAGCATCTCCCCGATGGCTCTGTGACTGTGGAGTCCATTCTAATCCAGGCTGCCGCCA
TTCGGAGGACCCTGGTACCAAGGTGCTATTGGTCTCCTGGACTTACCAGGATGAAGAGCTGGGAAGCTTC
CTCACAACATTGCTCAAAAAATGGCCTTCTCAGGCTCCCAGCTGACATGTATCCACCAGGAATCCACTGC
TCTTGTCTTTGCTGTTCTCCAGCAGGACCAGGGCTTCGGCCGAGAGGGGGTGCCTGAGTAAATGGG
AGGCCATGACTGGATTACAGACCCCGTTCTAGGGATAAGATCAGAGTCATAGGCTGGCTAAAGTAATAAA
CAGACCTTTTTGTTCTCTAAAAAAAAAAAAAAAA

Restriction Sites:	RsrII-NotI
ACCN:	BC008240
Insert Size:	909 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC008240 , AAH08240
RefSeq Size:	1222 bp
RefSeq ORF:	909 bp
Locus ID:	71957
Cytogenetics:	4 87.66 cM
Gene Summary:	Catalytic component of the Integrator complex, a complex involved in the small nuclear RNAs (snRNA) U1 and U2 transcription and in their 3'-box-dependent processing. The Integrator complex is associated with the C-terminal domain (CTD) of RNA polymerase II largest subunit (POLR2A) and is recruited to the U1 and U2 snRNAs genes. Mediates the snRNAs 3' cleavage. Mediates recruitment of cytoplasmic dynein to the nuclear envelope, probably as component of the INT complex.[UniProtKB/Swiss-Prot Function]