

## Product datasheet for **MC200620**

### LOC622757 (BC006659) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	LOC622757 (BC006659) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	LOC622757
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for BC006659, the custom clone sequence may differ by one or more nucleotides

```
CAGCAAGAAGCGGGTAGGCTTGCCTGACCTTAACCAGCCTGAAGCTTTCCTAGAAATAGGCTGGCAC
AGCTTGGCGGGCAATCGGTGGCCGCAATCCTCAGAAACTCCTTCAAAAAGAAGCTTGTCTTAGAGGCGT
GTGCACGTGACCTTGGCGCTCAGGAGGTGGAGGCGGGAGGATCAAAGGAGTTCAAAGCGTTTCCTTCTG
GGAAGTGAGTGGAGGAGGAAGGTCAGCTGGGTGCTTCTCTGCCTCTCTGACCTAGCAGGCTTGACCCCC
AGCATCTGGCCCCGTGAGTCTTCATTGGTATCGAACGGCTGGGATTTTTTGGTTTTCAAACAACATGTTG
GCTCTCTTAACTTTTGGCAAGACTCCGAGGGACACTGAGAGATCACTCCTGCTGAGGACAACCTGAGAAG
CCATCAGATTTGGTCCCCTGCATCAAGCCAGAAGCAGCCACAGGGATGGCTGTGTTTCTCATGACCCCA
ACAAGAAGACCACCATCTCACTGAGGCCAAGGAGTCAAGCTTGGTGTGGCACTCAAGCGCATATTGA
GGGCATCATGAAGAGGCCTCCTGAGGAGAGCTGCTGTTCAAGGACAACCAGCTCCTTGATGAGACGAAAA
CTCTGGGGGAATGTGGCTTCACTAGTCAGGTAGCAAAGCCACAGGCCCCAGCCACAGTGGCCTGGCCTT
CAGAGCAGAAGATGCCTTTGAAGGCCTGCCTATCTACCCCAATTCTTGTCTCCAGAACTCCAGGTGTG
ATGAAGGATTCTGGAAATAGTGCTAAGGAACAAGTTGTACAGTGAGGGCCCCAGGCCACCTCCTAGAGA
TCCATTCTCCCCACCACCACCACCAAAAAAGGGAGTTCGGGGCTGGTGTGAGTGGCTCAGTGGGT
AAGAGCACCCGACTGCTCTCCGAAGGTCCAGAGTTCAAATCCCAGCAACCACATGGTGGCTCACAACCA
TCTGTAATGAGATCTGACTCCCTCTTCTGGAGTGTCTGAAGACAGCTACAATGTACTTCCATATAATAAA
TAAATAAATCTTT
```

Restriction Sites:	RsrII-NotI
ACCN:	BC006659
Insert Size:	943 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC006659</a></u>
<b>RefSeq Size:</b>	1063 bp
<b>RefSeq ORF:</b>	943 bp
<b>Locus ID:</b>	622757
<b>Cytogenetics:</b>	16 A1