

Product datasheet for MC200718

Scamp4 (BC018215) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Scamp4
Synonyms:	2410022D05Rik; AI448940; AU022392; AW742369; Sc4
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC018215

```

CTTGGATAAGCGTGGCGTTTCGCTCGGGTGTCCCTCAGGTGGCCCCCTCTGCCTGGGTGCCAGGTCAGCTG
ATACAGAAATGGCAGGAAAGGAGAAACAACCTCCCCCACTGCCTCCCTTCTGCCACTGAAGCCCTGCTT
CTACCAAGGACTTCTCTGATGAGATCCCCGTGGAGCACCAGGTGCTGGTGAAGAGGATCTACCGGTGTGG
ATGTTTTACTGTGCCACTCTGGGTGTGAACCTGGTGGCTGCCTGGCTGGTGGATTGCGGGCGGAGCTG
GTGCCAATTTTGGCCTGGCCTTGCTCTGGCTGGTTCTCTTCACTCCCTGCAGTTATGTGTCTGGTTCCG
GCCTGCATACAAGGCCTTCCGAGCCGACAGCTCCTTTAACTTCATGACATTCTTCTTCATCTTTGGAGCC
CAGTTTGTGCTGACTGTATTAGGCCATTGGCTTCTCGGGCTGGGGCGCATGTGGTTGGCTGGCAGCTG
TTGGCTTCTTTGGGACCAAGTGTGGGGCTGCTGTGGTTATGCTGGTCCCGGCCATCCTGTTTTCTTGTG
CGCCCTTGTGATGGCCGTCACATATCGTGAAGGTGCACAGGATCTACCGTGGGGCTGGTGAAGCCTCCAG
AAGGCACAGACAGAGTGGAGCGCTGGTACCTGGAGGAACCCACCATCGAGGGAGGCCAGTTTAAACAGTT
TCTCTGGGAACAGCCTGCCTGAGTACCCCACTGTCCCAGCTACTCGTCCAGCGAGGCCACTGGCCTTA
GGCGGCTGCCTGGCCCCAGCACCATTGGCTCTGTCTCTGTTACCACTCGGTCTGAGGTCTGAGCATC
TCTCAAGGGTCCCAGGAGGTCTTGGCCCATGTTTCATCTCTTGCCCAAGCCCTCAGCAAGCACTGCCAA
CAACCATTTGCTCACTGCCATCCTACCAGAGCCTGGCTGGATTGCCCTGCTCCACCTGTTAAGAGGGTA
GAAACCTCAGCCTGTGCAAGGAGAGGCCTCAACCTCCAAGCGCCTGTCTTTGTGTTGCTGAGATGTAT
AAAGCCAACTCTGAGCGCTTGTGTAGATGGGTGTGAGTGGACAGTCAGGAGCATGGTGGCTGTGCCTGA
ATCCTCTGTTCTGCCATCCACTATACCTGTGCTCCTCTGGACCCTGCTACTTCTTGTCCCATCCCTTGT
GACAGTGACAACCTGCTGTCCCTCCACCACCCAGGATCCTGAATAGGCCCTTGGCCACTTGCAGACCT
GATCCACTTGGTTGAGGGTGTCTGGCTGCCGGGAGACCCCTCCCTCTCTGTGGGCAACTGGTGGGAGA
GGAACAGGCTCTGTACCTGGCAGGCCATACCCCAAGCCACAGGACTGCTGGGACAGATCTGTGCCTTAG
AGAGACCTGTACTTCTCCCAACCCCTCTGCCTTCCCGGGTAGCTGTGGCCTTGCCAGAAGAGCCAATC
TGAGTTTTGGAGTGTCCACAGCACCTCTGATGACATTGGACTCAATATTCTTCTTGTGACCCGATGTTA
CCAGCCAAGTCACCAAGGCTCCCCACAGGCTGTGTGCCAGTATCATGCCAACTGTTGGGGCACCTGCC
CTTACCCCATCGTGTGATGGAGTGGAGTCAAGTGCCCGCATGGGTAAACCCCTACAGCCACCTGCGG
CTCTTCCCTGTTCAACCCTGGGTATGTGCACCCACAGCCATAGTTGGTTTTGTATGTGAGAATATTCT
GGAAATAAATTCTCTTTCTGCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA AAAAAAA

```



Restriction Sites:	RsrII-NotI
ACCN:	BC018215
Insert Size:	693 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC018215 , AAH18215
RefSeq Size:	1828 bp
RefSeq ORF:	693 bp
Locus ID:	56214
Cytogenetics:	10 C1
Gene Summary:	Probably involved in membrane protein trafficking.[UniProtKB/Swiss-Prot Function]