

## Product datasheet for **MC200728**

### **Pglyrp1 (NM\_009402) Mouse Untagged Clone**

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Expression Plasmids  |
| Product Name:             | Pglyrp1 (NM_009402) Mouse Untagged Clone   |
| Tag:                      | Tag Free   |
| Symbol:                   | Pglyrp1  |
| Synonyms:                 | Pglyrp; PGRP; PGRP-S; Tag7; Tasg7; Tnfsf3l   |
| Mammalian Cell Selection: | Neomycin   |
| Vector:                   | PCMV6-Kan/Neo (PCMV6KN)  |
| E. coli Selection:        | Kanamycin (25 ug/mL)   |
| Fully Sequenced ORF:      | >BC005582 sequence for NM_009402<br>GGGCCAGCTCATACAGCCCTGCAGTCCTGTGCGGCACGTCCAGCATGTTGTTTGCCTGTGCTCTCCTT<br>GCCCTCCTGGGTCTGGCAACCTCCTGCAGTTTCATCGTGCCCGCAGTGAGTGGAGGGCCCTGCCATCCG<br>AGTGCTCTAGCCGCTGGGGCACCCAGTTCGCTACGTGGTGATCTCACACACAGCCGGCAGCTTCTGCAA<br>CAGCCCGGACTCCTGTGAACAGCAGGCCCGCAATGTGCAGCATTACCACAAGAATGAGCTGGGCTGGTGC<br>GATGTAGCTACAACCTCCTTATTGGAGAGGACGGTCATGTCTATGAAGGCCGAGGCTGGAACATCAAGG<br>GTGACCACACAGGGCCCATCTGGAATCCCATGTCTATTGGCATCACCTTCATGGGAACTTCATGGACCG<br>GGTACCCGCAAAGCGGGCCCTCCGTGCTGCCCTAAATCTTCTGGAATGTGGGGTGTCTCGGGGCTTCTTG<br>AGATCCAACATGAAGTCAAAGGACACCGGGATGTCAAAGCACTCTCTCCAGGTGACCAACTCTATC<br>AGGTCATCCAAAGCTGGGAACACTACCGAGAGTGAGAGACCTTGAGACCTAGTGAGAGTCCCCCCCACCC<br>CGAAATCCCTCCTGCCACCTGCTTCTTCCATTGACCCCAATAAAGACTCAGCACCAACAAAAAAAAA AAAAAAAAAA |
| Restriction Sites:        | RsrII-NotI   |
| ACCN:                     | NM_009402  |
| Insert Size:              | 549 bp   |
| OTI Disclaimer:           | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).   |
| Components:               | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).   |



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC005582](#), [AAH05582](#)

**RefSeq Size:** 713 bp

**RefSeq ORF:** 549 bp

**Locus ID:** 21946

**UniProt ID:** [O88593](#)

**Cytogenetics:** 7 A3

**Gene Summary:** Pattern receptor that binds to murein peptidoglycans (PGN) of Gram-positive bacteria. Has bactericidal activity towards Gram-positive bacteria. May kill Gram-positive bacteria by interfering with peptidoglycan biosynthesis. Binds also to Gram-negative bacteria. Involved in innate immunity. May function in intracellular killing of bacteria. The soluble form triggers apoptosis in vitro.[UniProtKB/Swiss-Prot Function]