

## Product datasheet for MC200729

### Higd1b (NM\_080846) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Higd1b (NM_080846) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Higd1b
Synonyms:	2310056K19Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC005568 sequence for NM_080846 AACAGATCCATGACCAAGAGGGCCCAAGATTCAAAGACTGTTCCAGTCTACTCTGGATTTTTCCAGCTG TTTTTCTTGGCCAGAGGATCACACCCATCGGGATGCAAACCGAGAGCAGGAAGCAAGCGGAACAGTAAAG TCAGATGAACCAGAGGAGAGGCAGATGTACGGACACGAGGACGGTAGCTTCGAGCAGCCCAAACCTTG CAAGAGAACACCCAGAGGGTATCTTTTCAAGTGCAGGTTGTGGTGGACATGGTCAGGACCTGAAGACTA AAAGTCAGAGTCTTATTGTGGCTATCCTCTTAGAAGTGGGCTAAAGCCTAGGCACACAGCTGTAA CATCTGGATCCAGGATCATGTCTGCTAACAAAGGGCTGGTGGTGCCACCTGAGGGTGAGGACAACCTGTC CAAGAAGTTCCTAAGGAAGACCAGGGAGTCTCCACTGGTGCCATAGGTGTAGCAGGCTGCCTGGTGATA GCAGCATACAGGATTTACCGGCTGAAGGCTCGGGGTTCTACCAAGCTCTCCATCCATCTGATTCACACCA GAGTGGCGGCACAGGCGTGTGCTGTGGGGCCATCATGCTTGGTGCCATGTATAAATGTACAGAGATTA CATCAAGAGGGTGTGAGGATGCCGAGAAGTAGGACTCCAGAGGTGGCAGGATGGCCAGTCCCCCGGC TGCTGCTGATGTTCTAATAAAGCTATGCTGAGGAAAAAAAAAAAAAAAAAAAAA
Restriction Sites:	RsrII-NotI
ACCN:	NM_080846
Insert Size:	297 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC005568](#), [AAH05568](#)

**RefSeq Size:** 753 bp

**RefSeq ORF:** 297 bp

**Locus ID:** 75689

**UniProt ID:** [Q99JY6](#)

**Cytogenetics:** 11 E1