

## Product datasheet for **MC200827**

### Slc50a1 (NM\_009057) Mouse Untagged Clone

#### Product data:

<b>Product Type:</b>	Expression Plasmids
<b>Product Name:</b>	Slc50a1 (NM_009057) Mouse Untagged Clone
<b>Tag:</b>	Tag Free
<b>Symbol:</b>	Slc50a1
<b>Synonyms:</b>	MmSWEET1; Rag1ap1; Rga
<b>Mammalian Cell Selection:</b>	Neomycin
<b>Vector:</b>	PCMV6-Kan/Neo (PCMV6KN)
<b>E. coli Selection:</b>	Kanamycin (25 ug/mL)
<b>Fully Sequenced ORF:</b>	>BC014292 sequence for NM_009057 CCACGCGTCCGGGCGGAGCTTCGCCTCCGCAGCTCCGGGTCTGAGCTACTGGAACCTCCGCGCTGGGCTA GGGAACCGGCTCCGGTGGGATGGAGGCGGGCGGCGTGGCAGACTCTTTCCTTTCTAGTGCCTGCGTGCTC TTCACCCTGGGCATGTTCTCCTACTGGCCTCTCGGACCTCAGGCATATGCAGAGGACACGGAGCGTGGACA ACATCCAGTTCCTGCCTTTTCTCACCACGGATGTCAACAACCTGAGCTGGCTGAGTTACGGAGTCTTGAA GGGAGATGGGACCCTTATCATCGTCAATAGCGTGGGGGCGTGGTTCAGACTCTTTATATCCTGGCATAT CTGCACTACAGTCTCAGAAGCATGGTGTGCTCCTGCAGACGGCAACCCTGCTGGCTGTCCTTCTCCTGG GTTATGGTACTTTTGGCTTCTGGTGCCGGACCTTGGGCGGCTTTCAGCAGCTAGGCCTCTTCTGTAG CGTCTTTACCATCAGCATGTACCTCTCCCACTGGCTGATTTGGCCAAGATCGTTCAGACTAAATCAACC CAGCGCTCTCCTTCTCCTGACCATTGCCACGCTCTTTTGTCCGCTCTTGGTCTATTTACGGGTTTC GCCTCCGAGACCCATACATCACGGTGCCCAACCTTCCAGGAATCCTCACCAGCTTAATCCGCTCCGGCT TTTCTGCAAGTACCCTCCAGAGCAAGACAGGAAGTACCGCTCCTGCAGACCTGACCCAGGCACCCGAG TGCCAACCTGGATACCAAACAGAGCTCCTGTTTCTGCTGGTCTTGTGACCAGTTCATGGATGCAATAGG TTGTGAGAAAAAGATGACTTTGAACTAAAGGGACCAAAGATCTTTCCTTAGATCAGATAGGACCTGTGG GATGAAATCACATTTTTACAGGAGATCATCTTTTCTCATTTTCGGAGGCTGAGGTGGTATTAGAATGTGCC TAAAAATAAACTGTTCCCAACCCCTGCCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA A
<b>Restriction Sites:</b>	RsrII-NotI
<b>ACCN:</b>	NM_009057
<b>Insert Size:</b>	666 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC014292</a></u> , <u><a href="#">AAH14292</a></u>
<b>RefSeq Size:</b>	1051 bp
<b>RefSeq ORF:</b>	666 bp
<b>Locus ID:</b>	19729
<b>UniProt ID:</b>	<u><a href="#">Q9C XK4</a></u>
<b>Cytogenetics:</b>	3 39.04 cM
<b>Gene Summary:</b>	Mediates sugar transport across membranes (By similarity). May regulate the expression of RAG1 a gene involved in V(D)J recombination.[UniProtKB/Swiss-Prot Function]