

Product datasheet for MC200868

Fth1 (NM_010239) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fth1 (NM_010239) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Fth1
Synonyms:	FHC; Fth; HfT; MFH
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	<p>>BC012314 sequence for NM_010239</p> <pre> CCGCGGTTTCCTGCTTCAACAGTGCTTGAACGGAACCCGGTGCTCGACCCCTCCGACCCCGCCGGCCGC TTCGAGCCTGAGCCCTTTGCAACTTCGTCTGTTCCGCCGCTCCAGCGTCGCCACCGCGCTCGCCCCGCCG CCACCATGACCACCGCGTCTCCCTCGCAAGTGCGCCAGAAGTACCACAGGACGCGGAGGCTGCCATCAA CCGCCAGATCAACCTGGAGTTGTATGCCTCCTACGTCTATCTGTCTATGTCTTGTATTTTGACCGAGAT GATGTGGCTCTGAAGAACTTTGCCAAATACTTTCTCCACCAATCTCATGAGGAGAGGGAGCATGCCGAGA AACTGATGAAGCTGCAGAACGAGCGAGGTGGCCGAATCTTCTGCAGGATATAAAGAAACGACCGTGA TGACTGGGAGAGCGGGCTGAATGCAATGGAGTGTGCACTGCACTTGAAAAGAGTGTGAATCAGTCACTA CTGGAAGTGCACAACTGGCTACTGACAAGAATGATCCCCACTTATGTGACTTCATTGAGACGTATTATC TGAGTGAACAGGTGAAATCCATTAAGAACTGGGTGACCACGTGACCAACTTACGCAAGATGGGTGCCCC TGAAGCTGGCATGGCAGAATATCTCTTTGACAAGCACACCCTGGGACACGGTGATGAGAGCTAAGCTGAC TTCCCCAAAGCCACGTGACTTTACTGGTCACTGAGGCAGTGCATGCATGTCAGGCTGCCTTCATCTTTTC TATAAGTTGCACCAAAACATCTGCTTAAGTTCTTTAATTTGTACCATTCTTCAAATAAAGAATTTTGGT ACCCAAAAAAAAAAAAAAAAA </pre>
Restriction Sites:	RsrII-NotI
ACCN:	NM_010239
Insert Size:	549 bp


[View online »](#)

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC012314 , AAH12314
RefSeq Size:	860 bp
RefSeq ORF:	549 bp
Locus ID:	14319
UniProt ID:	P09528
Cytogenetics:	19 6.23 cM
Gene Summary:	<p>Stores iron in a soluble, non-toxic, readily available form. Important for iron homeostasis. Has ferroxidase activity. Iron is taken up in the ferrous form and deposited as ferric hydroxides after oxidation. Also plays a role in delivery of iron to cells. Mediates iron uptake in capsule cells of the developing kidney.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) encodes the protein.</p>