

## Product datasheet for **MC200897**

### Nxph3 (NM\_130858) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Nxph3 (NM\_130858) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Nxph3  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >BC007167 sequence for NM\_130858

```
CGGAGGGGGCCTGGGACCCGAAAGGAAAAGAGAGACGGGTGAGAGCAAAGGAGGAAGATGCAACTGACC
CGTTGCTGCTTCGTGTTCTAGTGCAGGGCAGCCTCTATCTGGTCATCTGCGGCCAAGACGATGGCCCGC
CTGGCTCAGAGGACCCTGAACACGATGACCATGAGGGCCAGCCTCGGCCTAGGGTGCCTCGCAAGCGAGG
TCATATCTCACCTAAGTCCC GCCCTTGGCCAACCTCCACCTCCTAGGGCTGTTAGCCCCACCCGGGGAG
GTTTGGGGTGTCTTGGGCAGCCCCCAACCGCCCCAAACAAAGTCCCTTACCCTCGACCAAGGTGAAAA
AAATATTTGGATGGGGTGATTTCTACTCCAACATCAAGACAGTGGCCCTGAACCTGCTGGTCACAGGCCA
GATTGTGGACCACGAAACGGGACCTTCAGCGTCCACTTCGCCACAACGCCACGGGCCAAGGTAACATC
TCCATCAGCCTTGTACCCCCAGCAAAGCTGTGGAGTTTACCAAGAGCAGCAGATCTTCATCGAAGCCA
AGGCATCCAAAATTTTCAACTGCCGATGGAGTGGGAGAAGGTAGAACGGGGCCGCCGACCTCGCTCTG
TACCCACGACCCAGCCAAGATCTGTTCCCGTGACCACGCTCAGAGCTCGGCCACCTGGAGCTGCTCCCAG
CCCTTCAAAGTTGTCTGGTCTATATTGCCTTCTATAGCACGGACTATCGACTGGTGCAGAAGGTGTGCC
CAGATTACAACCTACCACAGTGATACCCCTACTACCCATCTGGGTGACCTAGATGGGCGGCAGAGGCCGG
TGCAGGAGGCTGGCTGGCTGGACTGGCAGGAAAAGGGATGGACCTCAGGTAGGGAGGGGGTGGAGAG
GAAGAGAGGCCTGGTGGGTGAGGATCAAGCGTCAGGTGTTGGGGAAGGGTCCCCAGTGCTGTGCCAACC
TGAGGTTAGAAGATACTTGTATGGAAGAGAAATATGCCCTGGCAGCCTTACAGGGCTTTCCCTCTAT
CCAAAGCAGATGCTGGGCCCGAGGCCCTGGGCAGGGAAGCCAGGCAACCCAGATCAAGCCACAGAGA
AGACCTCAACCTTTGTCTCTCTCTCTGCTATTGAGAGGAGGACAGTGGTCAGCAAAGGGGAGTCCAAC
GTGTAGACAACCTTAGGATGGAGAGCTCAATGGCGTCTCCCGTGGGAGCCGGCAGTGGGCAGCCCCTTCC
CCCCCCTCCCCCGCACTGAGGCTGAAATGGCAACCTCATCCCCACCTTGGGTGTTGCTGTCATGAC
TCTCGGTCAACAGCCAGGCCATTCTTTCCAGAAATCCCCTCTCCCTCCCTCCAGCACCTTTGTACCGCT
GACACCTCTTGCTGATGGCGCGCCACCCTTGTCTTCAGCTGAGACAGGATGACTATAGGTCTGTAGC
ATGAGGGGCACAAAGCAGAAGTCCGCCATGTGCTGTATGCCATCCTCTGCCCTAACCACTGTGAGA
GTGCTCCCTCCCTTGACAGCCAGTTCTTCTGTAAAGACTGGGGCATTCTGGGATAGGGATAGAGCCTCCT
TTGCTGGGCGGTGGGGCTATCCGCAGAGGAAGGCTCAGGAGGCTATGTGGCCATGCTGGAAGTGTGCCCT
AGTGTGTCTGTGGGGTGTGTGTGGGGGAGGGGTGGGAAGTCTGTGAGATGCCCTGTCACTTCTCTGT
GCAGAGTCTTGTCTTGGAGCAGGAAATAAAGCATTCCCAGGAAAAAAAAAAAAAAAAAAAAAAAAAAAA
```

**Restriction Sites:** RsrII-NotI



[View online »](#)

<b>ACCN:</b>	NM_130858
<b>Insert Size:</b>	759 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC007167</a></u> , <u><a href="#">AAH07167</a></u>
<b>RefSeq Size:</b>	1816 bp
<b>RefSeq ORF:</b>	759 bp
<b>Locus ID:</b>	104079
<b>UniProt ID:</b>	<u><a href="#">Q91VX5</a></u>
<b>Cytogenetics:</b>	11 D
<b>Gene Summary:</b>	May be signaling molecules that resemble neuropeptides. Ligand for alpha-neurexins. [UniProtKB/Swiss-Prot Function]