

Product datasheet for **MC201033**

Cyp2c50 (NM_134144) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cyp2c50 (NM_134144) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Cyp2c50
Synonyms:	AI266900
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC011222 sequence for NM_134144

```

CTCCATGGATCCAATCCTGGTCCTGGTGTTCACCCCTCTCCTGTCTGTTTCTCCTCTCACTCTGGAGACAG
AGCTCTGAAAGGGGAAGCTCCCTCCTGGCCCCACTCCTCTTCCAATTATTGGAAATATCCTTCAGATCA
ATGTGAAAGGACATCTGCCAATCCTTCACCAATTTATCAAAAGTTTATGGCCCTGTATATACTCTATATTT
GGGCAATAAGCCCACTGTGGTGTTCATGGATATGAGGCAGTGAAGGAAGCTCTGGTTGACCATGGGGAG
GAGTTTGTGGAAGAGGAAGGCTCCAGTTTTTGATAAAGCTACTAATGGAATGGGCATTATATTTAGTA
AAGGAAATGTATGGAAAAACACAAGGCGTTCTCACTCACGACTCTGCGGAATTTAGGCATGGGAAAAAG
GAGCATTGAGGACCGTGTTCAGAGGAAGCACGTTGCCTTGTGGAGGAACCTAGGAAAACCAATGGCTCA
CCCTGTGATCCCACCTTCATCCTGGGCTGTGCTCCCTGCAATGTCATCTGCTCAATTATTTCCAGGATC
GTTTTGATTATAAAGATCGGGATTTTCTTAATTTGATGGAAAACTGAATGAAATAACCAAGATAATGAG
CACCCCTTGGTTGCAGGTTTGAATACATTTCTGTTCTGCTTGATTATTGTCCAGGAAGTCATAACAAA
GTTTTTAAAAATTATGCTTGATTAAAAATTTCCTTTTGGAGAAAATAAAGAACATGAAGAATCATTGG
ATGTTACAATTCCTCGGGACTTTATTGATTATTTCTTAATTAATGGAGGCCAGGAAGACGGCAACTATCC
ATTGAAAAATAGGCTTGAACACCTGGCAATAACAGTGACTGATCTTTTTTCTGCTGGGACAGAGACAACA
AGCACAACTGAGATATGCTCTCCTACTCCTGTTGAAGTACCCACATGTACAGCTAAAGTCCAGGAAG
AAATTGAGCATGTGATCGGCAAGCACCGAAGCCCATGCATGCAGGACAGGAGCCGCATGCCCTATACTGA
TGCCATGATTTCATGAGGTCCAGAGATTCATTGACCTTGTCCTTCAACAGCCTGCCCCATGAAGTGACCTGT
GACATTAATTCAGGAATACTTCATCCCAAGGGGACAAATGTAATAACATCACTGTATCAGTGCTGC
GTGATAGCAAAGAGTTCCCAACCCAGAGAAGTTTGACCTGGGCACCTTTCTAGATGAGAACGGAAGTT
TAAGAAAAGTGACTACTTCATGCCCTTCTCAACAGGAAAACGGATATGTGCAGGAGAAGGCTGGCACGC
ATGGAGCTGTTCTATTCTAACCAGCATTTTACAGAACTTCAACCTGAAACCTCTGTTTACCCAAAGG
ACATAGATGTAACCCCAATGCTCATTGGATTGGCCTCAGTGCCTCCTGCTTCCAGCTCTGCTTCATTCC
TTCCTGAAAAATCAAATGCCTGGCTCCATCTGGGACATCTTCTGTGATCACCTGAGGCTTTAATTCAGTG
CTTTCACACTGGACATACTGCTTTTTATCCATTTCTTTTGCCTGCTCATAATTTATAAGGTTCTATGA
ACATACATTTACCAATTTGCAGATCTTGACAAATATTTCTGACTGGTCTATTTTCTGAAACACTGTTA
CCATATGCTGTGATGATGAGACCTGATTCCAAGTAAATATACTCTCATAAAATTGAAATTACTATTGT
AAGATAATTCAGTGTCGCTTCTCTCTGCATTTTCTTAATAATAAATTCTTTTGGTATACTAAAAAAA AAAAAAA

```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_134144
Insert Size:	1473 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC011222</u> , <u>AAH11222</u>
RefSeq Size:	1828 bp
RefSeq ORF:	1473 bp
Locus ID:	107141
UniProt ID:	<u>Q91X77</u>
Cytogenetics:	19 C3
Gene Summary:	Metabolizes arachidonic acid to several midchain and omega-terminal hydroxyeicosatetraenoic acids (HETE).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (1).