

Product datasheet for MC201065

Txn1 (NM_011660) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Txn1 (NM_011660) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Txn1
Synonyms:	ADF; AW550880; Trx1; Txn
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC010756 sequence for NM_011660 CCATTTCCATCTGGTTCTGCTGAGACGCGTGTGGCTCCCTCCCCGCAACAGCCAAAATGGTGAAGCTGAT CGAGAGCAAGGAAGCTTTTCAGGAGGCCCTGGCCGCCGCGGAGACAAGCTTGTCGTGGTGGACTTCTCT GCTACGTGGTGTGGACCTTGCAAAATGATCAAGCCCTTCTTCATTCCCTCTGTGACAAGTATTCCAATG TGGTGTTCCTTGAAGTGGATGTGGATGACTGCCAGGATGTTGCTGCAGACTGTGAAGTCAATGCATGCC GACCTTCCAGTTTATAAAAAGGGTCAAAAGGTGGGGAGTTCTCCGGTGCTAACAAGGAAAAGCTTGAA GCCTCTATTACTGAATATGCCTAATCATGCTCTGAAAAGTGTAACCAGCTACCAGCTGTTTAAACCTGT ACCTTTTTTAATTTGCAAAAACTATGAAGTGTGGAGAGTCTATACCAACTGCCATCTGATTATAAATG ACAATAAAATATTAATTCTACCCTTTTAAAAAAAAAAAAAAAAAAAA
Restriction Sites:	RsrII-NotI
ACCN:	NM_011660
Insert Size:	318 bp


[View online »](#)

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC010756 , AAH10756
RefSeq Size:	534 bp
RefSeq ORF:	318 bp
Locus ID:	22166
UniProt ID:	P10639
Cytogenetics:	4 31.87 cM
Gene Summary:	<p>Participates in various redox reactions through the reversible oxidation of its active center dithiol to a disulfide and catalyzes dithiol-disulfide exchange reactions (By similarity). Plays a role in the reversible S-nitrosylation of cysteine residues in target proteins, and thereby contributes to the response to intracellular nitric oxide. Nitrosylates the active site Cys of CASP3 in response to nitric oxide (NO), and thereby inhibits caspase-3 activity. Induces the FOS/JUN AP-1 DNA binding activity in ionizing radiation (IR) cells through its oxidation/reduction status and stimulates AP-1 transcriptional activity (By similarity).</p> <p>[UniProtKB/Swiss-Prot Function]</p>