

Product datasheet for MC201112

Farsa (BC013533) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Farsa (BC013533) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Farsa
Synonyms:	0610012A19Rik; Farsla
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC013533

```

AGCACGCTGGGAGTAGCCATGGCGGATAATCCCGTTTTGGAGCTGTTGCTGCGGCGGCTGGAGGTAGCTG
ATGGCGGTCTGGACAGCGCGGAGTTGGCTACGCAGCTGGGTGTGGAGCACCAGGCGGTGGTGGCGCAGT
GAAGAGCCTTCAAGCTCTGGGTGAGGTCAATTGAGGCTGAGCTTCTGTTCTACCAAGTGCTGGGAGCTGACT
ACTGAGGGTGAGGAGATTGCCCGGAGGGCAGCCATGAGGCCGTGTGTTTAGAAGCATCCCTCTCGAGG
GCCTGGTGAGAGTGAGCTCATGCACCTGCCAGCGGCAAGGTGGGCTTCCAGCAAGGCCATGTCCAACA
GTGGATCCGAGTGGACAAGAGTGCCGCCGACGGGCCCCGGGTGTTCCGAGTGGTGACAGTATAGAGGAT
GAAGTGCAGAAGCGTTGCAGCTGGTCCAGGCAGGCCAGGCTGAGAAGCTGGCTGAAAAGGAACGGAATG
AGCTCCGGAAGAGGAAGCTGCTGACGGAAGTGATCCTGAAAACCTACTGGGTGAGCAAGGGCAAGGCCTT
CAGCACAAGCGTGTCTAAGCAGGAGGCAGAGCTGAGCCCAGAGATGATCTCCAGTGGCTCTTGGAGGGAC
CGACCTTTCAAACCTACAACCTTCTCTGCCCGTGGTGTGCTCCCAGATAGTGCCACCTGCACCCCTTGC
TCAAGGTCGGCTCTCAGTTCGACAAATCTTCTGGAGATGGGGTCACTGAGATGCCTACGGACAATTT
CATTGAGAGCTCCTTCTGGAACCTTGTGCACTCTTTCAGCCTCAGCAGCACCCAGCACGTGACCAACAT
GATACTTTCTTCTGCGAGACCCAGCTGAGGCCCTGCAGCTTCCAATGGACTATGTCCAGCGCGTGAAAC
GGACCCACTCCCAGGGTGGCTACGGCTCACAGGGCTACAAGTATACCTGGAAGCTGGAAGAGGCCCGGAA
AAACCTGCTTCGCACACACACCACAGCAGCCAGCGCCCGGCCCTCTACCAGTTAGCCCAGAAAGCCGTTT
ACACCAGCAAGTACTTCTATTGATCGTGTGTTTCCGAATGAGACGCTGGACGCCACGCACCTGGCGG
GGAGTTCTTCAACAAGTGGGAATCACGCAGTTGCGCTTCAAACCAGCCTACAATCCCTACACAGAGCC
AGCATGGAAGTGTTCAGCTACCACCAAGTCTAAAGAAGTGGGTTGAAGTTGGCAACTCTGGGGTCTTCC
GCCCTGAGATGCTCCTGCCTATGGGGCTCCCTGAGAATGTGTCTGTCTATTGCTTGGGGTCTCTCCCTGGA
ACGCCCAACAATGATCAAATATGGCATCAATAACATCCGAGAAGCTGGTGGCCACAAGGTGAACCTGCAG
ATGGTATACGACAGCCCTGTGTGCCGACTGGATATTGAGCCAGGTCTCAAAGACACAGGAGGCTGCAT
GATATAGGCTGTCCAGAGATCCCTGTCTAGTCCCTGTGCATCCCTTGGCAGTGACCTACTATTTATGAG
GCCTCTGTGAGGCCAGCTCACTCCCCTGCCTCTCCCACCCTGGGGCAGGGTTCTAAGTATGGGGAGTT
CTTATGCTGGCTTTTGTGGTGTGTGTAGTTTGTGTAAGGGTGGGTTTGGCCCTGCAGGCCAGCTGTT
TTAATAAAGTGGGCACATTCCCTCGTTTGGTGCCTTTTGTAGGAAAAAAAAAAAAAAAAA

```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	BC013533
Insert Size:	1524 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC013533</u> , <u>AAH13533</u>
RefSeq Size:	1809 bp
RefSeq ORF:	1524 bp
Locus ID:	66590
Cytogenetics:	8 C2-C3