

Product datasheet for **MC201166**

Tmem176a (NM_025326) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem176a (NM_025326) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Tmem176a
Synonyms:	0610011I04Rik; AU040201; AU041743; Keg2
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC010831 sequence for NM_025326

```

GTGAGCTCTGGCCTTCCCTGCTGCACACTTCCCGGTCCCAGAGGAATCAGACGATGTCCACAGACATGG
AGACTGCAGTCGTCGGCAAGGTGGACCTGAGGCTCCACAACCCACCCACATTGATGTGCATCCACCA
GGAGTCTGCTCTGGCCAACTTCTGCTGGCCGGATGCTCATTGCTAAGGATTCCAGCATCCGCTTCCACC
CAGAGCCAGGGCAGCAGCAGAGTCTGGTGGCCTCCTGGGTGGTGCAGACTGTGCTGGGGCTCTGAGTG
TGGTTCTGGGTGGAACCTCTACATAGGCCATTATTTAGCCATGTATTCCGAAGGCGCCCTTCTGGAC
TGGGATCGTGGCTATGCTGGCTGGAGCTGTTGCCTTCTTCAAGAAACGGGGTGGTACCTGCTGGGCC
CTGATGAGGACCCTTCTGTGCTGGCAAGTTTCTGCACCGCTGTGGCTGCCATCGTTATTGGGTCTCGTG
AGTTGAATTTTACTGGTATTTTCTCGGAGATGATGTCTGTCAAAGAGACTCTTCATATGGATGGTCCAC
CATGCCTAGAACCACTCCAGTTCCTGAAGAAGCTGATAGGATTGCCTTGTGCATATACTACACAAGCATG
CTAAAGACCCTGCTCATGAGCCTCCAAGCTATGCTCTTGGGTATCTGGGTGCTGCTGCTCCTGGCTTCTC
TCACCCCTGTATGTGTCTACATCTGGAAGATTTTTCACAAAGGCGAAACAGAGGAGAAGAACTGCT
GGGTGCAGCTGTGATCTAGCCTTCTCTTCTGCTCCGGGCGTCCCTCCTACTGAAGCTGAAAGAAGAATC
AGGCAGGACTAAGAAGACCCTCCCCACTAGCAGGGCCATGGCCACTGCCTGGTTCTGCCAGCACCACA
GCAGCTCTCAGCAGCACTTGCTTGTCTCTCCATCCTTACCCTCCTATATCCCTCCTCAGGCAGCAACTT
GATAATAAACTCTCCTGTTATTGCTAAAAA

```

Restriction Sites:	RsrII-NotI
ACCN:	NM_025326
Insert Size:	735 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC010831</u> , <u>AAH10831</u>
RefSeq Size:	1024 bp
RefSeq ORF:	735 bp
Locus ID:	66058
UniProt ID:	<u>Q9DCS1</u>
Cytogenetics:	6 B2.3