

Product datasheet for **MC201274**

Arl4d (NM_025404) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Arl4d (NM_025404) Mouse Untagged Clone
Tag: Tag Free
Symbol: Arl4d
Synonyms: 1110036H21Rik; Arf4l; Arfl4; Arl5; ARL6; AW456149
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC016113 sequence for NM_025404
 GGCTGAGACGCGACGAGTCGGTGTCTGGTTGGCGTGGCTCGCGGAGAGTGC GGCTGAAGACACCCC
 AGCTTGTGCTACCTGCCGTTTTAAGTCTCCTCTTCTTGGCCTCTGCTGGAGTCCCAGACCTACTCCTGA
 GGTCCACCCTGCACCCTAGCTCACCATGGGGAACCACTTGACTGAAATGGCCCTACAGCCTCGTCTTT
 CTGCCGCACTTCCAGGCCTGCATGTCGTCTGATTGGTTGGATTCCGCGGGGAAGACTTCCCTTCTTT
 ACCGTCTCAAGTCAAGGAGTTTGTCCAGAGCGTCCCCACCAAGGGCTTAACTGAGAAGATCCGCGT
 GCCCTGGGGGTTCCCGTGGGATCACTTCCAAGTGTGGGATGTTGGGGTCCAGGAGAAGCTTCGGCCA
 CTGTGGGCTCTACACCCGCGGACAGATGGACTGGTGTGTTGGTGGACTCTGCAGAGACTGAGAGGT
 TAGAGGAAGCCAGAATGGAGTTACACCGGATCAGCAAGGCTTCCGACAACCAGGGGTGCCAGTGTGGT
 GCTGGCCAACAAGCAGGACCAACCGGGGCACTGAGTGCAGCCGAGGTGGAGAAGAGGCTGGCGTCCGA
 GAGCTGGCAGCTGCCACGCTACCCATGTGCAGGGCTGCAGTGTGTTGGACGGGCTGGCCTGCAGCCAG
 GCCTGGAACATTTGTATGAAATGATCCTGAAGAGGAAGAAGGCACCTCGGTGAGCAAGAAGAGACGGTG
 ACCCGAAGCTGCCCTACCTACCTGCTAGGGGTCAGTACACATTATAGTCAGGCAAGGCGCCGTGACCT
 CTCAGCCAGCCTGGGTGCTGGATCTGTCCACCGCAAAGGACTGAAGGAGCATAGGTGGCCCTGTGCGCCG
 CAGGGGGGAGGGAAGTGAAGTGTCTTCACTTGCCTTACCTGAAGAGGGTAGGGAGGGGCTGCAAGACC
 GGAGAAAGATGTAATGCTAACTCTACCTCCTACCTCGGCCCTGGCTCTCCTGTCTTCCCTGGCATT
 ATTTGGAGATATTTTTATCTTTCAGGGGAGGGTGTGTAGAAGGTTATCCAGAGAGTGCATTTGGGG
 GTGAATTGAGAACATCTCACAGTCGAGGAGAGGGTCTCCTCAGCTGCTTCTCACAGCTGTTCTGTG
 TCTGTGAAAATGCCAAGAATTCTCCCAACAATTGTAGATCCCCGGCCCTTTTATAAGCCTCGAGTGTCC
 TCTGTATTACTGTATTTTGTGTTTTTGGGTTTTTTTTTTTTTATCATTGCTGTAAGTTATTA
 AGCTGGATGATTCTGGTTAAAAAAAAAAAAAAAA

Restriction Sites: RsrII-NotI
ACCN: NM_025404
Insert Size: 606 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC016113 , AAH16113
RefSeq Size:	1363 bp
RefSeq ORF:	606 bp
Locus ID:	80981
UniProt ID:	Q99PE9
Cytogenetics:	11 D
Gene Summary:	Small GTP-binding protein which cycles between an inactive GDP-bound and an active GTP-bound form, and the rate of cycling is regulated by guanine nucleotide exchange factors (GEF) and GTPase-activating proteins (GAP). GTP-binding protein that does not act as an allosteric activator of the cholera toxin catalytic subunit. Recruits CYTH1, CYTH2, CYTH3 and CYTH4 to the plasma membrane in GDP-bound form (By similarity).[UniProtKB/Swiss-Prot Function]