

## Product datasheet for **MC201298**

### Krtap4-16 (NM\_001013823) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Krtap4-16 (NM\_001013823) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Krtap4-16  
**Synonyms:** OTTMUSG00000002196  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)

**Fully Sequenced ORF:** >BC016249 sequence for NM\_001013823  
GGAAATCTACTCAGAACCTCCATCTCTGACATCATGGTCAACTCCTGTTGTGGCTCTGTCTGCTCTGAGC  
AGGGCTGTGATCAAAGTCCCTGCCAGGAGAGCTGCGGCCAGCCAGCTGCAGTCAGCCAGCTGCAGCCA  
GCCAAGCTGTGCCAGACTACCTGCTGCAGGACTACCTGCTGCCGCCAGCTGCTGTGTGTCCAGCTGC  
TGCAAGGCCACCTGCCAGACCACCTGCTGCCGCCAGTCTGCTGCCAGACCACCTGCCGCCAGCTGTG  
GAGTGTCCAGCTGCTGCCGCCAGTCTGCTGCCAGACCACCTGCCGCCAGCTGTGGAGTGCCAGCTG  
CTGCCGCCAGTCTGCTGCCAGACCACCTGCCGCCAGCTGTGGAGTGCCAGCTGCTGCCGCCAGTC  
TGCTGCCAAACCACCTGCCGCCAGCTGTGGAGTGCCAGCTGTTGCCGCCAGTCTGCTGCCAGACCA  
CCTGCCGCCAGCTGTGGAGTGCCAGCTGCTGCCGCCAGTCTGCTGCCAGACCACCTGCCGCCAG  
CTGCGGAGTGCCAGCTGCTGCCGCCAGTCTGCTGTCAGACCACCTGCTGCCGTACAACCTGCTGTGG  
CCCTCTTGCTAAATTTCTCTAACTACCAATCCCGAGGCAGCCATAATCACAATGCTCTAACCTTCTC  
ATCCATGTAACTTTTTCTGTACTAAATGCCTCCTACCTTTGTTGCTCCTGACACCAGTAAGAACTGCC  
AAAAGAGAGTGCTACATTACGTACCGCTCATATTTCCCAAAGATCTCACCTCTACATGTCTGTCTCAGT  
ACATGGTATAGCTGATCATATCTCTCAAGTCTCTACTACCATCCTTCCCTTGGCCTTTACCATTTCTGTG  
TAACTATTCTTCAATAAATATTGCCTTATTTTCAGAAAGTATAAAAAAAAAAAAAAAAAAAAAAAAAA

**Restriction Sites:** RsrII-NotI

**ACCN:** NM\_001013823

**Insert Size:** 609 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC016249</a></u> , <u><a href="#">AAH16249</a></u>
<b>RefSeq Size:</b>	972 bp
<b>RefSeq ORF:</b>	609 bp
<b>Locus ID:</b>	435285
<b>Cytogenetics:</b>	11