

Product datasheet for MC201317

Lman1l (NM_199222) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lman1l (NM_199222) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Lman1l
Synonyms:	BC020188; slamp
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC020188 sequence for NM_199222

```

CCCACGCGTCCGCCACGCGTCCGGTCTCCCGCCTGACTGCGTGGCGTCTATAGCCCATCGTTGGCCT
CACAGCCAGGATCCGTCCAGGTGCCAGAAACCAGCTCCTTTGCAATGCTGGAGATCAGAGGCCTCTCAC
CTTCCCTCTGCCTCCTTTCCCTGTTGCTGGTTCTGCACGGTGCTGAAAGGAGCCAGCCCCCTCCTCGACG
CAGGTTTGAGTACAAGCTCAGCTTCAAGGGACCAAGACTGGCAGTTCTTGGGGCTGGGATACCCCTCTGG
AGCCACCATGGAGACGCCATCCTAGGCTTAGAGGAAGTGGCGTGGTTCCATCCATGAAGAACGGGAGCG
GTGCCGTGTGGAGCAACATCTCTGTCTCCTTCCCTTCTGGGAGGTAGAGATGCAGATGAGAGTGACTGG
ACCAGGACGCCGGGAGCCAGGGTGTGGCCATGTGGTACACCAAGGACAGGGCTCAAGTCGGCTCTGTT
GTAGAGGAGCTAGCCTCATGGGATGGCATCGGGATCTACTTTGATTCTCCACCAGTGATGTCCAGGACA
GCCCTGTCTCCAGTGCTGGCCAGTGATGGGCATGACCTCCAGGAGCAGTCTGGGGATGGAATGTCCG
GGAGCTTGGCTCCTGTCTACAGGACTTCAGGAACCGGCCATTCCCTTCCGAGCTCGGGTACCTACTGG
AGGCAGAGGCTGCGTGTGTCTTTGAGTGGAGGCCCTTACTCCTAAGGACCTGAAGAGGTCTGTGTTGATG
TGGAGCCCTGTTTTAGCCCTGGAGGCTTCTTTGGGGTCTCAGCAGCCACTGGCACCTTAGCAGCAGA
TGACCATGACGTCTGTCTTCTGACCTTCAGTCTGAGGGAACCAAGTCCAGAGGAAACCCCCAGCCA
TTCATGGAGAAGGATCAGTTGCTCCTGGCCAGGAAGCTGGAAGAACTGAAGGCAAGGCTGGCCCTGGGCA
CCAGAGAGGCTAGCATACCACCACTGAATCCCAAAGCCAGGAAGAAGGTGAAAGGTTCTTTAACCTGGA
GGATACGCTGGGAGACAAAGCCAGATCCTACAGGCCCTCCAGGCTCTCTCCAGGAGATGGCCAGGCT
GAGAAGCACTGGAACAGCAGCTGGGGTCTACAGTCCAAGTCCGGCCTGAGGGAGGCTGGGAAGCAGCT
GCCCGGATGGCTTCAGGAGCACAGCTTCTACCTGCCTGTGGGCACTAAGCATCACTTTTTTGAAGTGG
ACCATATCCTCAGCCTCCTGCAGAAGGACCTCCGGGGTCTGGTGCAGGCAGGAGCTGGACAAGAGGCTTC
AGGAGTACCTGTCCACAGGAAGCCTTTCTAGAGCCTGCATTACCCATCACCAGGACGATAGGGTTCT
GAGGAGACAGCCATCTCCCCAGCATGAAGCCTGACCCAGCTCAAAGCCATGTCTTACTGGGAAAAA
AAATGGGGATAGAGATTTAGCCAGTGCCACATCTCAGTTTTTGGCAAGTCCATCATATTAAGGTGA
TTTTCTCTTGCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAAAAAAAA

```

Restriction Sites: RsrII-NotI



[View online »](#)

ACCN:	NM_199222
Insert Size:	1125 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC020188</u> , <u>AAH20188</u>
RefSeq Size:	1629 bp
RefSeq ORF:	1125 bp
Locus ID:	235416
UniProt ID:	<u>Q8VCD3</u>
Cytogenetics:	9 B