

## Product datasheet for **MC201517**

### Mrps6 (NM\_080456) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Mrps6 (NM\_080456) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Mrps6  
**Synonyms:** AW046321  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)

**Fully Sequenced ORF:** >BC027271 sequence for NM\_080456  
GGAGCTGACCGGCGTGTTCCTCTCGTCTTCTCCTCTGCGCCCCGCTCCCCGCCCTCGCGACCCGGCT  
CTCCTGGACTCGGCGCCCAACCTGGGCGATGCCCGCTACGAGTTGGCTTTGATTCTGAAAGCCATGC  
GGCGGCCAGAGACCGCTGCTGCTTTGAAACGTACAATAGAATCCCTGATGGACCGAGGAGCCATAGTGAG  
GAACTTGAAAGCCTGGGTGAGCGTGCGCTCCCCACAGGATCTCGAGTCACAGCCAGCAGCACAGCCGA  
GGAGGGTATTTCTGGTGGATTTTATGCTCCGACAAGTCTGTGGAGAACAATACTGGAACACTGGCGC  
GAGACATTGACGTGGTTAGACCAAATATTGTGAAACACCCTCTGACCCAGGAAGTAAAAGAGTGTGACGG  
CATAGTCCAGTCCCACTGAAGAAAACTGTATTCAACAAGAGGAGGAAGAAGTGAAGAAGATTACCA  
GATTCTGGCCTTATTTAATCCTAAGGGCACTATGGGTGCTGCTAGGTTGTTGCTAGGATACTTTAGC  
CCATGACCATTTGCTGCAGGAGGTAGAACTGCTGGCCGAGACCTGCCCTGATGTCTCTGCTGAGATTT  
CATCCCACTTGTGGGTTTGTGCGGAGTGGGGTGTTCACAGTACCACTGTAGCGTTTCCAAGAGCAAAA  
TGTTTGTCACTTCACTTGGTTGCTTGAAGCCTATATGGAACACTGGGAGCAGAGTAATAAATCATGAC  
TTTATCAACAAAAA

**Restriction Sites:** RsrII-NotI

**ACCN:** NM\_080456

**Insert Size:** 378 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online >](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC027271](#), [AAH27271](#)

**RefSeq Size:** 795 bp

**RefSeq ORF:** 378 bp

**Locus ID:** 121022

**UniProt ID:** [P58064](#)

**Cytogenetics:** 16 C4