

## Product datasheet for MC201520

### Gng11 (NM\_025331) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Gng11 (NM_025331) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Gng11
Synonyms:	0610037B21Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC039062 sequence for NM_025331 ACCTGGTCCTGTTCTCCTGCAGGGTCTGCAGCCTGCCAGAGCAGCCAGGAGGATGTCCCATATTTAGGCAG ACTGCCCTGGCTACAGGGCTGCTGCGGCCCTGGAGAAGCCGCCCTGAAGTTGGTTGGAGTCTGGGAAAGGCC ATCGGGAGCAGGCTTGGGCAGGAAGCCGAGCCCCGGGGAGCTGGGAGCTGCGGTTCCCGAGCCAGTCAC GTCCTGCGCGCGGGGGCGGTAGCCTGGGCCCGCGGGTATTTGGGGAAGCGTCCCAGTTCTAGGCGCTTC CAGTGGCTCCCTGCACAAGTGACTTCTGGGCCAGTTCTGGGGCAAAGATGCCTGCCCTTCACATCGAG GATCTGCCGAAAAGGAAAACTGAAGATGGAGGTTGAGCAACTTCGAAAAGAAGTCAAGTTGCAGAGAC AACAGGTATCTAAATGTTCTGAGGAAATAAAGAACTACATTGAAGAAGTCTGGAGAGGATCCTCTGGT AAAGGGAATCCCAGAAGACAAGAATCCCTTCAAAGAAAAGGGCAGCTGTGTCATTTATAAGTAACCTCG GGGGAGGGAAATCTGTTCTTTTGTAAATTTTCCCAAATGAAGCCAAAGTGTGTGTTGAAGAACTGA GAAAAAATAAAATCGAGAAAAGACTGTCATATAAGCACTCTCCAAGCAGTTGGGGTGTGAATCATACT GCTTCTCACCTTGGCTCCCCACACTCTCTCAGAGGAGAGCATCTGATTGCAGTTATGGAATGAAGG CATGGCTCCCCACTTTTCACTACTCCTGTCTCCAATAAAGTTTTATCTTGGAAAAATCAAAAAA AAAAAAAA
Restriction Sites:	RsrII-NotI
ACCN:	NM_025331
Insert Size:	222 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC039062](#), [AAH39062](#)

**RefSeq Size:** 849 bp

**RefSeq ORF:** 222 bp

**Locus ID:** 66066

**UniProt ID:** [P61953](#)

**Cytogenetics:** 6 A1

**Gene Summary:** Guanine nucleotide-binding proteins (G proteins) are involved as a modulator or transducer in various transmembrane signaling systems. The beta and gamma chains are required for the GTPase activity, for replacement of GDP by GTP, and for G protein-effector interaction (By similarity).[UniProtKB/Swiss-Prot Function]