

Product datasheet for **MC201597**

Fam167b (NM_182783) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fam167b (NM_182783) Mouse Untagged Clone
Tag: Tag Free
Symbol: Fam167b
Synonyms: SEC
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC030183 sequence for NM_182783
 CCACGCGTCCGCATCTGCTCACTCAGCCTGGCACACTTAGCGTCCCTCACCCACCTTCTCAATTCAGCAC
 ACACCTTTGGCACATTCTTTTCAACCTCATCTTTTCACTTCTCACACGGCCAAATAGCCAGCCTTTCC
 CTTACTGTCCCGTTGGTCACCACAGACCTACCATGTCCCTAGGGCCGCTGAAATTCAGGCAGTGGGAGA
 GGAGGGTGAGGAGGATGAGGAAGAGAGCTTAGATTCTCTGAAGGCCCTGACTGCCAAGCTTCAACTGCAG
 ACCCGACGGCCCTCATACCTCGAGTGGACTGCCAGAGTCCAGAGCCGGGCTTGGTACAGAGCCCAAGCTA
 GACCTGAGCCAGTCGGACCCGGGCCATCTGTGGCTTTGACTCCATGGATTCTGCGCTTGAGTGGCTCAG
 ACGGGAGCTGCAGGAGATGCGGGCTCAGGACCGGCAGCTGGCAGGGCAGCTGCTCAGGCTGCGGGCCCGG
 CTGCACAGACTGAAAGTGGACCAAGTCTGTACCTGCACCAGGAGCTTCTGGATGAGGCTGAGCTGGAGA
 TGGAGTTAGAGTCTGGGACTGGCTTGCCTCTGGCCCCACCGCTGCGGCATCTGGGACTCACGCGCATGAA
 CATCAGTGCCAGACGCTTACCCTCTGCTGACAGCAGACTTGGGTGTCTCTTGCAGTATCTGGGGAGAAA
 GAAGGAAGGAAGAGGGACCCCGGAGGCTCTGGCTACCTGCTGGGAAGGTGGGCACACTTAGGTTTCCAA
 AAGCTGAATTTAGAGAGCACAGGATGGAGGGGAGGAGAGGAAACTCGGGCCCCAAATGTCTTAATAA
 AAAATGCATTGAATCCCAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: RsrII-NotI

ACCN: NM_182783

Insert Size: 225 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC030183](#), [AAH30183](#)

RefSeq Size: 874 bp

RefSeq ORF: 225 bp

Locus ID: 230766

UniProt ID: [P17257](#)

Cytogenetics: 4 D2.2