

Product datasheet for **MC201653**

Mgst3 (NM_025569) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mgst3 (NM_025569) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Mgst3
Synonyms:	2010012L10Rik; 2010306B17Rik; 2700004G04Rik; AA516734; GST-III
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC029669 sequence for NM_025569 CGGACCGTGGGCGAAGGTGAGCCAGCCAAGATGGCTGTCCTCTCTAAGGAGTATGGATTTGTGCTTC TCACTGGTGCTGCCAGCTTTGTGATGGTGCTCCACCTAGCCATCAACGTGGGCAAAGCCGCAAGAAGTA CAAGGTAGAGTACCCTGTATGTACAGCACAGATCCTGAGAACGGGCATATGTTCAACTGCATTCAGCGC GCCACCAGAACACGTTGGAGGTGTACCCTCCCTTCTGTTTTTCTAACGGTGGGAGGTGTTACCACC CGCGCATAGCTTCTGGCCTGGGCTGGCCTGGATTATTGGGCGAGTCTTTACGCATATGGCTACTACAC AGGAGACCCTAGCAAGCGGTATCGAGGAGCCGTGGGCTCTCTTGGCCTCTTTGCCCTGATGGGCACCACC GTGTGCTCTGCTTTCCAGCATCTCGGCTGGATCAGACCAGGCTTAGGCTACGGTCCAGATCCTGCCACC ACTGAGGTGTGGAGGGCCTTCCGACTCTACTCACCTCCAGCGACTCACCTGATTTCCAGTTGCACTGG TTTTTTTTTTTTTTAATAATAAAAACTTATCTGGCAAAAAAAAAAAAAAAAA
Restriction Sites:	RsrII-NotI
ACCN:	NM_025569
Insert Size:	462 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC029669](#), [AAH29669](#)

RefSeq Size: 614 bp

RefSeq ORF: 462 bp

Locus ID: 66447

UniProt ID: [Q9CPU4](#)

Cytogenetics: 1 H2.3

Gene Summary: Catalyzes oxydation of hydroxy-fatty acids. May participate to the lipid metabolism.
[UniProtKB/Swiss-Prot Function]