

## Product datasheet for **MC201718**

### Tspan31 (NM\_025982) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Tspan31 (NM_025982) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Tspan31
Synonyms:	2700085A14Rik; AW742554; Sas
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC023306 sequence for NM_025982 GACGGTCCTCGATACCCTCCCCACAAGCTTTGGGACCATTGGGTCGCCAGATCTGGGGAGATGGTTTGC GGCGGATTTGCCGTGCCAGGAATGCGCTGTGCGCTCTCAACGTGGTTTATATGCTTGTGGGCTTTTTGC TCATTGGCGTGGCTGCCTGGGGCAAGGGCTCGGTGTGGTGTCTAGCATTACATCATTGGAGGAGTCAT TGCTGTGGGAGTCTTCTCTCCTCATCGCGGTGGCAGGACTGGTGGGCGCCGCAACCACCATCAAGTG CTGCTCTTCTTTTATATGATCATCCTAGGTTTGGTCTTTCATCTTCCAGTTTGAATCTCTTGCTCATGTC TGGCTATTAACAGAAACACACAGGCAGATGTCATCAATGCTTCGTGGTGGTCTTGAGCAACAGTACCCG GCATGAAGTGGAGAGAAGTTTTGACTGCTGTGGCTTGTCAACCTTACAACCCTGCGCCTACAAGATGAT ACTTCCTGTAGTGGGTGTGCAAACTAAGAGCTCTACATGCCAGATGTGTGGGAGAGGTTCTTAAAC ATTCAGACAAAGCCCTCAAGATCCTAGGAGGTGTTGGGCTCTTCTTTAGTTTACAGAGATCCTTGGTGT TTGGCTAGCAATGAGATCCGGAATCAAAGGACCCTCGAGCCAACCCAGTGCCTTCTATGAGACTTT TGACCTCCAGATAGCATCTCTGCAACCTAGGCTGTTTCTAGCTAGCTCTGTATCCAAGGAAGACCTTGA ATTCCAGACCCTCCTGCCTCAGCCTCCTGCTGGCATCACCACTCAGCCAGCCCTGGATCCTTCTGAC TTGGCTTTGCCATTTTTTTGTTTTTTGAGACAGGGTTTCTCTGTATGGCCCTGGCTGTCCTGGAAGTGTG TAAACTAGGCTGGCCTCAAACCTCAGAAATCCACCTGCCTCTGCCTCCCAAGTGTGGGATTAAGGGCGTA CACCACACTGCCTGGCCTTGCATCTTACAACACCTAGGGGTGGCCCTATGGTTTTGTTTTAGAGAA AAGAGTCCCTTACCTCACCCCAAGAAGTGGTGGTGGACTTATGTCTGTTTTGATGGGAACACGATTA TGGATGAGGGGAGGTAGGCCTCGGGCCCTAGGAGGTAAGTGCACCCCAAGACCCTTGGCTTCA AGGACAGAGTCAGAACCAAGAGAAAACCAAGTTATAACCCAGGGCCTCACAAGTGTCTGGCAAGCATT CTACCAATCAGTTGTCTTCCCAGCCACCTACAGTCTCATCAAAGCAGGACCAGTGAACCTTAGAAACT AAAGCCATTTTGACCAATATTCAAAGAAAGCATTCTCTCTGATCTTGCCAGCCAAAGAGCTGGACTGC CTGCCTAGGCTCTGCAAGGAAGTGTGTTTTGTTGTTTTTTTTGTTGTTGTTTTGTTTTCCCAAGT CTCCAGTGGCAAAGCCAATTGGAAAAACATTAAGAAAAATAAGCCCAAAAAAAAAAAAAAAAAAAAAA AAAAAAA
Restriction Sites:	RsrII-NotI
ACCN:	NM_025982



[View online »](#)

<b>Insert Size:</b>	633 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC023306</a></u> , <u><a href="#">AAH23306</a></u>
<b>RefSeq Size:</b>	1548 bp
<b>RefSeq ORF:</b>	633 bp
<b>Locus ID:</b>	67125
<b>UniProt ID:</b>	<u><a href="#">Q9CQ88</a></u>
<b>Cytogenetics:</b>	10 D3