

Product datasheet for **MC201760**

Lad1 (NM_133664) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lad1 (NM_133664) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Lad1
Synonyms:	lad-1; LADA
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC025580](#), [AAH25580](#)

RefSeq Size: 2558 bp

RefSeq ORF: 1587 bp

Locus ID: 16763

UniProt ID: [P57016](#)

Cytogenetics: 1 E4

Gene Summary: Anchoring filament protein which is a component of the basement membrane zone.
[UniProtKB/Swiss-Prot Function]