

Product datasheet for MC201960

Folr2 (NM_008035) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Folr2 (NM_008035) Mouse Untagged Clone
Tag: Tag Free
Symbol: Folr2
Synonyms: FBP; FBP2; Fol; Folb; Folbp-2; Folbp2; FR-; FR-beta; FR-P3
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC022108 sequence for NM_008035
AAAGCTTCTATGTGGGGCTGTGGACGAAGACTGTAGAGACTACCCAGAGTCTGACCTAGGGAGAGGCCAA
CTCGGATACCCCTATGTGCGCTCCCAGAAGCTAAGGACATTGAGACAGAAAGACATGGCCTGGAACAGA
CACCACTCTTGCTTTTGGTCTACATGGTACAACAGGCAGTGGCCGGGACAGAACAGACCTACTCAACGT
TTGCATGGATGCCAAACACCATAAGACAAGCCGGGCCCGAGGACAAGCTGCATGACCAGTGTAGTCCA
TGAAGAAAAATGCCTGTTGCTCAGTCAACACCAGCCAGGAGCTACACAAGGCTGACTCCCGTCTGACT
TCAACTGGGATCACTGTGGCAAGATGGAGCCTGCCTGTAAGAGTCACTTCATCCAAGACTCCTGCCTGTA
TGAGTGTCCCCAACCTTGGGCCTTGGATCCAGCAAGTGGACCAGAGTTGGCGTAAAGAGCGTTTCCTG
GATGTGCCCTTATGCAAAGAGGACTGTCACCAGTGGTGGGAAGCCTGTCTACCTCTTTACCTGCAAGA
GAGACTGGCATAAAGGCTGGGACTGGTCCCTCAGGCATTAACAAGTGCCCAAACACAGCACCTGTACAC
GTTTGAGTACTACTTCCCGACACCAGCCAGCCTTTGCGAGGGTCTCTGGAGTCACTCCTACAAGGTCAGC
AACTACAGCAGAGGGAGTGGCCGCTGCATCCAGATGTGGTTTGACTCAACCCAGGCAATCCAATGAGG
ACGTGGTGAAGTTTTATGCTTCCCTTTATGACATCTGGGACTGTGCCCCATGCAGCAGTACTTCTGTGCC
CAGCCTGGCCCCAGTGTGTCATTATGGCTCCCTGGCTGAGAGGTCAGTCTTCTCTAGATTTCTCT
CTATCTACCCCTGGTCTGGTTCAACTCTTCAAAGAATAAGGAAGTCTTGAGCCTGCTCCACCCCTCTCC
TCTGTCATCCAGTTCTGATCCATGTTGGGGTTGGGGTTTCTACAATCATTTTCAATAAATCTATGACA
CATCTGGAAAAAAAAAAAAAAAA

Restriction Sites: RsrII-NotI
ACCN: NM_008035
Insert Size: 756 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC022108</u> , <u>AAH22108</u>
RefSeq Size:	1072 bp
RefSeq ORF:	756 bp
Locus ID:	14276
UniProt ID:	<u>Q05685</u>
Cytogenetics:	7 E2
Gene Summary:	<p>This gene encodes a receptor protein located on the plasma membrane that mediates folate uptake by cells. Mice lacking the product of this gene show no defects in embryonic development and grow normally into fertile adults. However, such mice were found to be highly susceptible to the teratogenic effects of arsenic. Alternate splicing of this gene results in multiple transcript variants. [provided by RefSeq, Dec 2014]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. All variants (1-3) encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>