

Product datasheet for MC201972

Slc35a3 (NM_144902) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc35a3 (NM_144902) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Slc35a3
Synonyms:	2310050P13Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC024110 sequence for NM_144902

```

CCCACGCGTCCGCCACGCGTCCGCCCGAAACGCCTGGGCCTCACCGAAGCGGTCCGGGCCACAGAGCC
GGCGGGAGCCGGCGGTGCTGCGGGTCAGCGGACGCGGCGAGCCGACGGGCTCGGCAGCCTGGGTCTCCC
GGCGTCGCGGAGCCCGGGGCTCCAGTGCGGGCAGCTGGCCGAGCCGCACGGCGGACAAATGAAGAGAAAA
CAATGTCTGCCAACCTAAATATCTTTCTTGGGAATTTTGGTGTTCAGACTACCAAGTCTGGTTCTAAC
GATGCGGTATTCTAGGACTTTAAAGAGGAGGGCCTCGTTATCTGTCTTCTACAGCAGTGGTTGTGGCT
GAATTTTTGAAGATAATGGCCTGCATCTTTTAGTCTACAAAGACAGTAAGTGTAGTGTGAGAGCACTGA
ATAGAGTACTGCATGATGAAATCTTAATAAGCCCATGGAACCCCTGAAGCTCGTATCCCGTCAGGAT
ATATACTTTTCAAGAACTTACTCTATGTGGCACTGTCAAACCTAGATGCAGCCACTTACCAGGTTACA
TATCAGTTGAAAATACTTACAACAGCATTATTTTCTGTGTCTATGCTTGGTAAAAAATTAGGTGTGTACC
AGTGGCTCTCCCTAGTAATTCTGATGGCAGGAGTTGCTTTTGTACAGTGGCCTTCAGATTCTCAAGAGCT
GAACTCTAAGGACCTTTCAACAGGCTCACAGTTTGTAGGCCTCATGGCAGTTCTCACAGCCTGTTTTTCA
AGTGGCTTTGCTGGAGTTTATTTTGAGAAAATCTTAAAGAAACAAAACAGTCAGTATGGATAAGGAACA
TTCAACTTGTTTTCTTTGGAAGTATATTTGGATTAATGGGTGTATACGTTTATGATGGAGAATTGGTCTC
AAAGAATGGATTTTTTTCAGGGATATAATCAACTGACGTGGATAGTTGTTGCTCTGCAGGCACTTGGAGGC
CTTGTAATAGCTGCTGTCATCAAATATGCAGATAACATTTTAAAGGATTTGCGACCTCCTTATCCATAA
TATTGTCAACAATAATATCTTATTTTGGTTGCAAGATTTTGTGCCAACCAAGTGTCTTTTCTTGGAGC
CATCCTTGTAATAGCAGTACTTTCTTGTATGGTTACGATCCCAACCTGCAGGAAATCCCACTAAAGCA
TAGTCATGTCTGCTCCTTAACCTGATTCTCACAGTGGTGCCCTAGGAACCTCAACATTAATCTTGCACAG
AGGACTTCTGCAGATTCTATGAAGAAATCATCAGGCTAAATCTGATCGTTCTGTTCAAATGGTTTTGAAA
CGTGTAGCCTTCAGGATGAAGTGTGTCCGTCATCTGCAACCCAAGCTTGAAGAGTGTCCAGGAAGTAGG
CTTAGAGGGGATACTGGGAGGAGCCTTACGTCTGAGCTCCTAGTGACTTCAGCCTTTTGTGTGTAG
CACTGTCCAGTGCAGCCTTTGCAAGAGCACACATGCAAGAAAGAGCAAGGTTTGCAAGGTAGATCTG
CCTACACTTATCCCATGTTCTTAACCTTACCCTTTCTGTACAGATGTATCAAGCTACATGATCTATTTAA
ATGGAAATTTATAATTGCCTATCAAGCTGTGAACAAACAAATATAATTTAAACCAAAAAAAAAAAAAAA
  
```

Restriction Sites: RsrII-NotI


[View online »](#)

ACCN:	NM_144902
Insert Size:	981 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC024110</u> , <u>AAH24110</u>
RefSeq Size:	1680 bp
RefSeq ORF:	981 bp
Locus ID:	229782
UniProt ID:	<u>Q8R1T4</u>
Cytogenetics:	3 G1
Gene Summary:	Uridine diphosphate-N-acetylglucosamine (UDP-GlcNAc) transporter in the Golgi apparatus. May supply UDP-GlcNAc as substrate for Golgi-resident glycosyltransferases that generate branching of diantennary oligosaccharides (By similarity).[UniProtKB/Swiss-Prot Function]