

Product datasheet for MC202001

Emc4 (NM_026519) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Emc4 (NM_026519) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Emc4
Synonyms:	2610318K02Rik; Tmem85
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC051926 sequence for NM_026519

```

CTTGCAAGTTGCAGCTAGTTCTTACCTGCAGAACGCAAGCGGCCGAGGGTGTGGGTACTTTTAGCCAGCG
GCCGCTGTCATGACGACCCAGGGGGGCTAGTGGCCAACCGAGGCCGGCGTTCAAATGGGCCATTGAGC
TGAGTGGGCCGGGAGGAGGCAGCAGGGGCCGAAGTGACCGGGGCAGTGGACAGGGAGACTCCCTCTATCC
AGTCGGTTACTTGGACAAGCAAGTGCCTGATACAGCGTGCAAGAGACCGACCGGATCCTGGTGGAGAAG
CGCTGCTGGGACATCGCCTTGGGTCCCCTCAAGCAGATTCCCATGAACCTCTTCATCATGTACATGGCAG
GCAACACCATTTCCATCTTCCCTACTATGATGGTGTGTATGATGGCCTGGCGACCCATTCAAGCACTTAT
GGCCATTTTCAGCCACTTTCAAGATGCTGGAGAGTTCAAGTCAGAAGTTTCTTCAAGGCTTGGTCTATCTC
ATTGGGAACCTGATGGGTTTGGCATTGGCAGTTTACAAGTGCCAGTCCATGGGACTGTTGCCTACACATG
CATCAGACTGGCTAGCCTTTATTGAGCCCCCTGAGAGAATGGAGTTCAGTGGAGGAGGCTTACTTCTGTG
ACCCCCAGAAAGCACCATGCATATGGACCTCATCTGCATTTTAAGCTGTTGACTCCTAAATATCCTCCCT
TAAACAGTCATGCATGTGGAGAAGAGCCTCCTGTGGTGAGGCCAGTGTGCAGCATGTTACAGCAGAAAA
GAAACTACAGCACAACCTATGACTGAAAAATAATGTAGAAAACTTTATTTATGTTCCAGCAGAGCAAAA
ACAAACAAAACCAGACTCTATGTAAACAAAAGAAATAGCAGCTGCTCAGTAGCAGCTTCTCCTTTCAATAA
ATTAAACGGTTGAGAACAATGAAAAAAAAAAAAAAAAAAAAAAAAA
  
```

Restriction Sites:	Ascl-NotI
ACCN:	NM_026519
Insert Size:	552 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC051926</u> , <u>AAH51926</u>
RefSeq Size:	955 bp
RefSeq ORF:	552 bp
Locus ID:	68032
UniProt ID:	<u>Q9CZX9</u>
Cytogenetics:	2 E3
Gene Summary:	May mediate anti-apoptotic activity.[UniProtKB/Swiss-Prot Function]