

Product datasheet for **MC202024**

Prl7a1 (NM_008930) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Prl7a1
Synonyms:	AA407021; AA409971; AA410116; PLP-E; PLP-G; Prlpe; Prlpg
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC051671 sequence for NM_008930 <pre> AAGACAAGTGTGCAAGGACACCTCAGAGATGCCGCTGTCTTTCACTCAACCATGCTCCTGGGCATTCTG CTGCTGGTGGTGTCAAACCTCCTTTTATGGGAGAATGTGGCCTGTCTACCTTTAAGCAGCAATGATACTG ATGATGATCCACTATCCATCAAGGACTGTTGGATCATGCCATGATACTTTCTAAGAGTATCACTGACCT CAACATGGAGTTGCGCAGGATATTTACCATCAGTGAGATGTCAGCTAAACTTATTGATAAAATTTCTAAGT TCATCATCATCATCAGACTCCTATGATCAATTTATGCTTGAATTTCTTGGGCAGCAGGAGTTACTGACTA AGAACCTCAGCTACTGCCACAAATACTCCATCAAAGTTCCCGAAGACATAGAAGAAGCCCAAAATGTCAT CTCTCTTGAAGACTTTCCAATCTTGATACTCAGTAGAATGCAGGCTTGGAAATGAAACTCTGAAAAACCGA ATCAACCTATCTGAAGGTACACCAGGAATAGACGATGATATCCTCCCTATATATAAAATATTGAGACAA AAATTGCAGAACTTCTTGAGGACAGCAAGAGTATACTCAGCCAGGCTTATGGAGCAACAGAAAATGTGGC TGATTACACCTCTGGTCTGGTCTTGAAGACTTTCAATCATCTGATGAAGAACTCGATTTTGGCTCTT TGTAATTATCCTATTGTTTGCATGTTGATATCCACACAGCTAACTTTTATCTCCAGTTCTTGAGGTGTG TGGCTCTTGTAAATAGTGACAGCTGCTTATCTTCCAAACTGGAATGATTGATGCTGTATTATTTT AAAATAGTCTGATTTTATGCATTTCAAAGATGAGCATGAGTAAATGGGCATCTTTTAAAGAAAAATAA ACATTGTGTGTTTAAACAAAGTTTAAAAATTAATAAAAAAAAAAAAAAAAAAAAAA </pre>
Restriction Sites:	RsrII-NotI
ACCN:	NM_008930
Insert Size:	798 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC051671</u> , <u>AAH51671</u>
RefSeq Size:	959 bp
RefSeq ORF:	798 bp
Locus ID:	19113
UniProt ID:	<u>O54830</u>
Cytogenetics:	13 12.56 cM