

Product datasheet for MC202148

Ubxn8 (NM_178648) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ubxn8 (NM_178648) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Ubxn8
Synonyms:	D0H8S2298E; Rep-8; Rep8h; Ubxd6
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC024492 sequence for NM_178648

```

CGGACGCGTGGGTCGCCCACGCGTCCGCGCGCCATGGCTTCTCGTGGCGTCGTGGGCCTTTTCTTGCTCT
CTGCTCTCCCTCTCTTGTGTCTGGAGCTCCGGCGCGGGATCCCGAGTCTCGGCATTAAAGACCTCATTTT
GCTGAGTGGCCGGATATTCTACTGCTGGCCCTGCTTACCTTAGTCATTTCCGTGACCACCTCCTGGTTT
AACTCACTGAAACCTTCGCAAGGTCATCTGAAAGAAGGAGAAAAAGAGAATGAGAAAAAGACGGAGACTTG
TGAGGGAAAGGCAACAGGAAGCCAGGGTGAGAAGGCCAGCAGATATATAGAGAATGTTTTAAAGCCTCA
GCAGGAAATGAAGTTGAAAAAAGTGGAGGAGCGGTTTTATCAGATGACAGGGGAGACCTGGAAGCTGACC
GCTGGCCACCGACTCCTAGAGGGGGATGAAGATTCGGAGTTTGAAAAATCAAGTCAGGCATCTTTTGAAA
CAATAAAGCAGGAGGCAGCAAGAAGACAGAACTTGCCTAAGTTTTCAACTGAAATTTACCAGCTGCTCG
GCCGCTCCTGAGGAAGGAGGTTCTGATTTACCTGAAGAACCATCAGAAACAGCTGAAGAGGTAGTTACT
GTTGCCCTTCGATGCCCCAATGGACGTGTTCTGAGGAGAAGATTTTTCAAGTCGTGGAATTCACAGGTCT
TACTTGACTGGATGATGAAAGTTGGGTACCACAAATCCCTATATAGACTTCCACTTCCTTTCCAGAAAG
GGCCCTGGAAGTGAGGGGAGGCTCGTCACTGGAGGACATAGGAATAACTGTGGATACTGTACTTAACGTG
GAAGAGAAAGAGCAGAGCAGCCAGTAGAGAGGCTGGAAGCACCTGCTGCCCTTCCACAGGGCGGTGTTCA
TTCAAGCCCCGATGCACCTGGAGCTCCTTGTCCACAGGTGTACCCCGCCACACCCCGCCACACCCC
ATGGATGAGGACCGCTTTGCACTGATGGTCTGTGGCCTGTATTATAATGTGAAAGGAGGTAAAAAGAGT
AAGACAGACCTGGATGCAGGCCTCAGCCGGCATGGTGGGACACGCCTCAACTCCCACAATTTCTGAGTT
AATTAGACCATCTTCAAACCAAAACCAACCCTAGCTCATATAAGAACTGGTTCCTGCTTGCGTCTACCTT
ACCAAAAGTATTTCACAATTCAGTCAGAATCTAAGGAGGATGAAGGGACATTCTATTGTATTTCTAAAA
ACTTAAAGAACTCCAGAGCCAACTGATCCTGTATAGTTTTGACTTAGATCTGGGGATGGAGGTCTGGTC
CATGAACCTTGGGGAACATAGTGAAACCACAGTAAACAAGGCCATTGTTGTTAAGTGCATTTGGAATGTG
AGCTGTCGTTTGTCTACAAATTGATGCCCTTAGAGCATCTCTCAACCTGTGGGGTCTGTGGGGCTCAAA
TAACCTGTTTATAGGGGTGCACCTAAGACCATCTGCATGAGATTTACATTGATGATTATAACAGCAGCA
AAATTGCTTATAAAGTAGCAATAAAATAATTTTATGGTTGGGGGGTCAACACAATATGAGGCAATGTATT
AAAGGGTTGCAGTATTACAGAGGTTGAAAACCTAAAAA
  
```

Restriction Sites: RsrII-NotI


[View online »](#)

ACCN:	NM_178648
Insert Size:	834 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC024492</u> , <u>AAH24492</u>
RefSeq Size:	1666 bp
RefSeq ORF:	834 bp
Locus ID:	108159
Cytogenetics:	8 A4