

Product datasheet for **MC202221**

Rmdn1 (NM_025476) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rmdn1 (NM_025476) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Rmdn1
Synonyms:	Fam82b; RMD-1
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC025158 sequence for NM_025476

```

CGGACGCGTGGGCTGTCAGTCATGGCGCTGTCTCGCGACTGTGGCGTCGTCTCCCTCCACTCCGCGTCT
GTCAGCAGACGAGGACCTGGAGCAGCCGTGGGCGCTGTGGCGCCTGGGGAGTCCGCGCTTCCGAGGTGAT
AGGAAATACAAGACCTTTCAAGAGAGGCTTTTTATTCTCAGCCTTGCTTATTTGGGTTTTGAAACCTAT
CAAATCATTTCTCAGGCTGCTGTGGTTCATGCCACAGCCAAAGTTGAAGAAATACTTGCGCAAGCAGACT
ACCTCTATGAAAGTGAGAAACAGAAAAGCTTTATCAGTTGCTTATACAGTACAAGGAAAGTGAGGATGG
AGAGCTACTATGGCGGCTGGCAGGACATCCCGTGATATTGCACAGCTTAGTAAGACTTCAGAAGAAGAA
AAGAAAGTTTTGGTGTATGAAGCTCTAGATTATGCAAAAAGAGCACTAGAGAAAAGGAGTCAAGCTCTG
CTGCTCATAAGTGGTATGCAATCTGCATTAGTGATGTTGGAGATTATGAAGGCATCAAGGTTAAGATTGC
AAATGCCTATGTTATCAAGGAGCATTTTGAGAAAGCAATTGAGCTGAACCCTAAAGATGCTACTTCAATT
CACCTTATGGGTATTTGGTGCTATACATTTGCTGAAATGCCTTGGTATCAAAGAAGAATTGCTAAAGTGC
TGTTTGCAAAATCCTCCTAGTTCCACCTATGAGGAGGCCTTAAGATACTTTCAAAAGCAGAAGAAGTGGA
TCCAAACTTCTACAGCAAAAACCTTGCTGCTTTTAGGAAAGACATACTTGAAGCTAAACAACAAAAAGCTT
GCTGCGTTCTGGTTGGTAAAAGCCAAGGGTTATCCAGCACACAGAGAAGATAAACAGATACAGACAG
AAGCTGCTCAGTTGCTTACAGGTTTGTGACAAGAAGTGAGAACATTTTGGAGAAGGCAACAACATACCTA
ATACTTTGTTTTCACTGATTTATAAACATTATACATTAACAAAAA
  
```

Restriction Sites:	RsrII-NotI
ACCN:	NM_025476
Insert Size:	918 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC025158</u> , <u>AAH25158</u>
RefSeq Size:	1031 bp
RefSeq ORF:	918 bp
Locus ID:	66302
UniProt ID:	<u>Q9DCV4</u>
Cytogenetics:	4 A3