

## Product datasheet for **MC202362**

### Polr2j (NM\_011293) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Polr2j (NM_011293) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Polr2j
Synonyms:	14.5kDa; Polr2i; Rpb11a; Rpo2-4
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)

**Fully Sequenced ORF:** >BC042939 sequence for NM\_011293  
GCGAGAGCGCGGCTAGGATGAACGCTCCTCCGGCCTTCGAGTCGTTCTTGCTCTTCGAGGGCGAGAAG  
AAAATCACCATTAACAAGGACACTAAGGTTCCCAACGCCTGCTTGTTCACCATCAACAAAGAAGACCACA  
CTCTGGGGAACATCATTAAATCGCAGCTGTTGAAGGACCCGAGGTGCTGTTTGCTGGCTACAAAGTCCC  
TACCCCTTGGAGCACAAGATCATCATTGTCGAGACACCCAGACTACAGTCCCCAGGAGGCTTTC  
ACCAACGCCATCAGACCTCATCAGCGAGCTCTCCCTTCTGGAGGAGCGATTCCGGGTGGCCATCAAGG  
ACAAGCAAGAAGGAATTGAGTAGCAGCTGAAAGAAGCGTTGCTTAGTGGCTGGAAGGCTGGCACATACTC  
CTCAGGGCCCTTCAGTTTACCACATGGCGACAGCCACTCTCAGCTCTTAGCTGGGCACGGGGCTTGCC  
CCATTCTGTACATAGACTGTTGTTTCTAATAAAAAGCAGAGCAGAAAAAAAAAAAAAAAAA

**Restriction Sites:** RsrII-NotI

**ACCN:** NM\_011293

**Insert Size:** 354 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC042939](#), [AAH42939](#)

**RefSeq Size:** 550 bp

**RefSeq ORF:** 354 bp

**Locus ID:** 20022

**Cytogenetics:** 5 G2