

Product datasheet for MC202691

Plin2 (NM_007408) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Plin2 (NM_007408) Mouse Untagged Clone
Tag: Tag Free
Symbol: Plin2
Synonyms: AA407157; Adfp; ADPH; Adrp
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC054766 sequence for NM_007408
 CTTGCTCGTCCCTCAGCTCTCCTGTTAGCGTCTCTTTTCTCCAGGAGGAAAAAATGGCAGCAGCAGTAG
 TGGATCCGCAACAGAGCGTGGTATGAGAGTGGCCAACCTGCCCTTGGTGGAGCTCTACCTACGACCTTGT
 GTCTCCGCTTATGTCAGTACAAAGGATCAGTACCCGTATTTGAGATCCGTGTGTGAGATGGCCGAGAAG
 GGCGTGAAGACCGTGACCTCTGCGGCCATGACAAGTGCCCTGCCATCATCCAGAAGCTGGAGCCACAAA
 TTGCGGTTGCCAATACCTATGCCTGCAAGGGGCTAGACAGGATGGAGGAAAGACTGCCTATTCTGAACCA
 GCCAACGTCGAGATTGTTGCCAGTGCCAGAGGTGCCGTAACGGGGCGAAGGATGTTGGTACGACTACC
 ATGGCTGGAGCCAAGGATTCTGTAGCCAGCACAGTCTCAGGGTGGTGGATAAGACCAAAGGAGCAGTGA
 CTGGCAGCGTGGAAAGGACCAAGTCTGTGGTCAATGGCAGCATCAATACAGTTTTGGGGATGGTGCAGTT
 CATGAACAGTGGAGTAGATAATGCCATCACCAGTCGGAGCTGCTGGTAGACCAGTACTTCCCTCTCACT
 CAGGAGGAGCTGGAGATGGAAGCAAAAAGGTGGAAGGATTTGATATGGTTCAGAAGCCGAGCAACTATG
 AACGGCTGGAGTCCCTGTCTACCAAGCTCTGCTCTCGGGCTTATACCAGGCTCTCAGCAGGGTTAAGA
 GGCCAAACAAAAGAGCCAGGAGACCATTCTCAGCTCCACTCCACTGTCCACCTGATTGAATTCGCCAGG
 AAGAATATGCACAGTGCCAACCAGAAAATTCAGGGTGCTCAGGATAAGCTCTATGTCTCGTGGTGGAGT
 GGAAGAGAAGCATCGGCTACGACGACACCGATGAGTCCCCTGTGTTGAGCACATCGAGTCACGACTCT
 GGCTATCGCCCGAACCTGACCCAGCAGCTCCAGACTACATGCCAGACTGTCTGGTCAACGCCCAAGGG
 TTACCACAGAACATTCAAGATCAGGCCAAACACTTGGGGTGATGGCAGGCGACATCTACTCCGTATTCC
 GAAGGAATCCTTAGATGAAGTTATGGATTACTTTGTTAAACAACACGCCTCAACTGGCTGGTAGGTCCC
 TTTTATCCTCAGTCTACCGAGGTGAACAAGGCCAGCCTGAAGGTCCAGCAGTCTGAGGTCAAAGCTCAGT
 AAACCCCTCCTTGTACCCAGAGCATGATGTTGCTGGCCAGATGACCCCTTTTGCTGTATTGAAATTAAC
 TGGTAGATGGCTTTAGCTTAGAAAAGCAGCTTCTTAGAACCAAGGGCCTCATTATGGTCACTCACAGCTC
 AGTTATGGTCTTGCCCCAGCTGGCCCTGGCACAGGAGTTCTTACCTGGCTGGTGGAGTGGCCTGTGTTA
 GTCTTGTGAGGACCTCGAGGAACCTAAAAGCTCAGATGCACCTACAGTCTTGTCTGTGGCCTTTGTATTG
 TTATTGGCTGTAACGTCTGTCTGGACCGAATAAAGATTCAATCAGTGAIAAAAAAAAAAAAAAAAAAAAAA AAAAAAAAAAAAAA

Restriction Sites: AscI-NotI



[View online >](#)

ACCN:	NM_007408
Insert Size:	1278 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC054766 , AAH54766
RefSeq Size:	1693 bp
RefSeq ORF:	1278 bp
Locus ID:	11520
UniProt ID:	P43883
Cytogenetics:	4 40.69 cM
Gene Summary:	May be involved in development and maintenance of adipose tissue.[UniProtKB/Swiss-Prot Function]