

Product datasheet for **MC202795**

Nudt10 (NM_001031664) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nudt10 (NM_001031664) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Nudt10
Synonyms:	AI853080; Dipp3a
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC055771 sequence for NM_001031664

```

GCCGGGCCCTCTCCCTGCCAGCCTTTGCTCCTCGGCTGCGCGGCGGCAGTGGCCCGCTCGCTAGGCAGCT
GCAAGCGACCTCAGAGGCGGCGGCGGCGGCGGCGTGGGGAGGATGAAGTGCAAGCCGAACCAGACGCGG
ACCTACGACCCCCGAGGGCTTCAAGAAGCGCGCCGCGTGCCTGTGCTTCCGACGAGCGCGAGGACGAGG
TGCTGCTGGTGAGCAGCAGCCGCTACCCCGACCGCTGGATCGTGCCGGGAGGGGGCATGGAGCCCGAGGA
GGAGCCGGACGGCGCGGCGGTGCGCGAGGTGTACGAGGAGGCGGGAGTCAAGGGGAAGTTGGGCCGCTTG
CTGGGCGTCTTCGAGCAGAACCAGGACCGCAAGCACCGGACCTACGTGTTCTGTGCTCACCGTACCGGAGC
TGCTGGAGGATTGGGAAGACTCGGTGAGCATCGGCAGGAAGCGCGAGTGGTTCAAGATCGAAGATGCCAT
CAAGGTCTGCAGTGCCATAAGCCCGTGACGCGGAGTACCTGGAGAACTGAAGCTGGGCGGCTCCCT
ACTAATGGAACTCGGCCGCCCGTCCCCGCCAGAGAGCGAGCCCTAGTGAATGGAATAAATGTTATTGA
GATTTACTTTGAGAGAATCAAGTATGTGAAGCCAGTGATGAATGGAATTGTGAAGCACAGGTGAGCTCTT
TGACCCCCCAAGGACACAGCTCATCCTATGCTTTTGGACACTTCTTTCCTGGTATTAACAATGGCTACT
CTTCAGGTTGTTGCACTACCACTGTATCTGTACAGAATTCTATTTGGCACCTGTGGCCTTTGTGTGCTTT
TGTTTTGATTTCCCATCAGATTTGACCGTGAATGCAAAGAGCGTGCTTTCTTAGTCTTTTATATCTAGC
ATACATTACTTGATATATGGTAGCGTCACTGAATATTTGTTTAAATTTTGAATGATAATGAATGGTG
AAACATCCCCCTTCTAATATCTCTAGCATATATGTGGCTTTTGTATTTCAAGGCAGCTGTGAGATAGGAG
CAAAAATGTTGTTTCATATGTCTGTGTTATCCAGACATTATGTATCTCCTACATTTGCCAAAATTTTAC
CTAATTTGACTGTGCCTCTTGACATAATTTCTCCTCTATAGAAAATGCAGGAGACAGAGAAACAGTTTA
CTATATATCAACAGATGAGATTCATGGAATCTTCAGACTCCAAGAAAACACAGGACAGATGATCCAGTT
TCTTCGGTATATTAGTCACAAGGATGAAAGGAAAAGGTGGAGGCTGAACCTGTGGGTTGAGAAAAGGCTTC
AGAGAATGATCAGTACCTATATGGATCTTATCTCGATCCTGATTTGATCAAACAGAAATGAGTGAATTT
ATCAGACAACAGAACTTTGAGCATTGGCTGGATATTTTCATGTTGATTATTATAAATACTATTGTTAT
ATTGGAATCAAGAATACTTATCTGCATTCAAATATTTACAGTTGAAATGATATAATGTCTGGGATTTGCT
TCAAAATATGTCTTGATATTAATACACTGGACATAACTCAAACAAGATTGGCCTTTAAGATCATCTAAG
CTCAGTGATAATTAACGATGTTTATTATATGTTTCTCTGCTTTTGTATGTTTGACATTTCCATA
ATAAAAAGTTAATCAAAATCATTGGAAAAA

```



[View online »](#)

Restriction Sites:	Ascl-NotI
ACCN:	NM_001031664
Insert Size:	495 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC055771</u> , <u>AAH55771</u>
RefSeq Size:	1729 bp
RefSeq ORF:	495 bp
Locus ID:	102954
UniProt ID:	<u>P0C027</u>
Cytogenetics:	X A1.1