

Product datasheet for MC202824

Lrrc42 (NM_029985) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lrrc42 (NM_029985) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Lrrc42
Synonyms:	A930011F22Rik; AA536996
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC027203 sequence for NM_029985 <pre> CCTCCCTGCTGGTACCTCAGTCCCGCGGCCAGGTGCGCCCCAGCCCCCTCCCTCTTGCCACCCTCTGTG TGTGTCCACCCCTCAGGGCTTTGTGGCTAGGTTGCCGAGCCGCTTCCCGGTCAGCCGCTCAGTGGAGG AAGGAGCCGCCAGTGAGGTGCAGCCGCTCGCTGCCCCGGGCAGGGGCACAGTTGCAACCAAGGCAATGGC TTACTACCTCAACTCAGAAGCCACCTGGATCCAGGGCCATCTATGTTGAGAAAAATGGGCAGCTACAC ATGGTCAACCTGGCTCTGGATGGTGTAAAAACAGCCTGCAGAAGCCACGACCATTAGACTCTTTCCCA AAGGCTTTTCCGTGGAGCTGTGCATGAACAGGGAAGACGACACTGCACAGAAAGAGAAGACCGATCACTT TATCTTCACGTACACCCGAGAGGGGAATCTGCGATACTCTGCTAAGTCCCTCTTCAGTCTGGTGTGGGC TTCATCTCCGACAATGTTGATCACATTGATTCCTCATCGGCTTTCCTGAGCAGATCGCTGAGAAGCTGT TCTCTGCTGCTGAAGCCAGACAGAAATTCACAGAGCCAGGGGCAGGGCTCAGGGCGTTACAGAAGTTCAC TGAGGCCTACGGGAGTCTGGTGTGTGCTCGCTGTGCTTGCGAAACAGGTATTTGGTGGTTGCAGAAAAG CTTGAAGAGATCAAGTCTTCCGGGAGCTGACCCGTCTAGACCTCTCCTGTTGCTGGCTTGAGATGAGC ATGAGCTTCTAGAACACCTCACCAATGAAGCTCTGTCCAGTGTCAGCTCCACTTGAAGGATAATTG CTTATCTGATGCAGGCATTCCGAAAAATGACAGCACCAGTGCGGGTGATGAAAAGAGGCCTTGAGAATCTA GCGCTGTTAGACCTGTCTTGTAACCCCTGAGATAACAGACGCTGGAATTGGATACCTGTTTTCATTAGAA AATTAAATTGCTTGATATCTCTGGAACAGGACTAAAGGACATCAAAGCCGTTAAGGACAAGCTCCGGGC AAACATAGGCCTTGTTCACTCCAAAGTGCCCTTTGAAGGAATTTGATCATAGTAAGTCAAAACAGAGGGC TGGGCTGACCAGATTGTTCTACAGTGGGAGCGAGTGAGCGTGGAAGCCGTGCGGCAGCGGAAAGATGTGG AGCCTCGAAAGGCAGCTCAATACTTCTACCAGAATCGAGCTCGACAGAAGTTACAAGGAAATGTCCCTC CGCAGAGACCCACATGAACCTCTTCTGGAAGAACTCCAGTTCTATAGGGAGGAGGCTCCGGACTGCCATGAG CCACTGTTGAGCCAAGAATCAAAGAAGAGCAAGAAGAGAGCTTTCAAGGAGTCAGAGCAGGAGCAGAGCT CGCCCCAGTCGGCCAAGCAGAAGTGTGTGTGCTGGCTGTGGAAGATTGGGATTTGCTCAATTCCTACTG ACTGACACATACAGTCGACCTCCCCACCTCAGTACGAGGAGGGGAGACTGCTGGTGATTTTCTTTTGGT TTGGATTTACCTATGGCATTAGCATAGCAAGCAGATATTTCTACTTTTGTGGTGCAGGAGGGAAGGCTAT GGACGTATGTAATAATTTCTAATGTATATTTTCAATACATAGTAATTTTTCAAATAAAAAAAAAAAAAA A </pre>
Restriction Sites:	RsrII-NotI


[View online »](#)

ACCN:	NM_029985
Insert Size:	1266 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC027203</u> , <u>AAH27203</u>
RefSeq Size:	1681 bp
RefSeq ORF:	1266 bp
Locus ID:	77809
UniProt ID:	<u>Q8R2U7</u>
Cytogenetics:	4 C7