

Product datasheet for **MC202990**

Mxra8 (NM_024263) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mxra8 (NM_024263) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Mxra8
Synonyms:	1200013A08Rik; 1700095D18Rik; AI131686; Asp3; Dicam
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >BC026438
 GCGGGTGTGGCAGGGAGAGGGGGTTCGGCTGTTTTCTGGAGAGTTAGAGCTGAGTAAGACAAAGCACGT
 CCCCCGACGGCGCCATGGAACCTGCTGTCCCGTGTCTGTGTGAAACTGCTGCTTCTCAGAGCTCTGC
 AGTCTGTCTCAGGGCCTTCAGGGACCGCAGCAGCCAGCAACTCTGTGTGTCTGAGTCTGTGGTGAGC
 TTGGCAGCCGGAACCCAGGCTGTGCTACGCTGCCAGAGCCCCGCATGGTGTGGACCCAAGACCGGCTGC
 ATGATCGCCAGCGCGTGGTCCACTGGGACCTCAGCGGGGGCCGGGCAGCAACGGCGCCGACTTGTGGA
 TATGTATTTCGGCGGGTGAACAGCGCGTGTACGAGCCGCGCGATCGCGACCGCCTCCTGCTGCGCCTTCT
 GCTTTCCACGACGCAACTTCTCGCTGCTCATTGCGCTGTGGAGAGAGCGATGAAGGGGTGTACACCT
 GCAACCTGCACCATCACTACTGCCACCTCGATGAGAGCCTGGCTGTGCGCCTCGAGGTTACAGAGGATCC
 CCTATTAAGTCGCGCATACTGGGACGGTGAAGAAGTGTGGTGGTGGCCATGGCGCGCCGGCACTG
 ATGACCTGCATCAACCGTGCACAGTGTGGACTGACCGCCATTTAGAGGAGGCGCAGCAGGTGGTCCATT
 GGGACCGACAGCTACCTGGAGTGTACACGACCGCGCGACCGCCTGCTTGACCTGTATGCATCTGGCGA
 GCGCCGCGCCTATGGCCGCCCTTCTGCGTGATCGCGTGTGAGTGAACACCAACGCTTTTGACGCGGT
 GACTTCTCCCTACGCATCGATGAGCTGGAGCGAGCTGATGAGGGCATCTATTCTGCCACCTGCACCATC
 ACTACTGTGGCCTCCACGAGCGCGAGTCTTCCACCTACAGGTACAGAGCCTGCCTTTGAGCCACCAGC
 TCGTGCTTCTCCTGGCAATGGGTCTGGTACAGAGTGTCTCTAGCCAGATCCACCCCTGACCCGAGGC
 CACAGCATCATCAATGTCTTGTCCCAGAGGACCACACATTTCTCCAGCAACTGGGGTACGTGTTGG
 CCACGCTGCTGCTTTCATCTTGTCTGCTCATCACTGTAGTCTGGCTACACGACATCGTCACAGCGGAGG
 ATGCAAGACGTCGGACAAAAAGCTGGGAAGTCAAAGGGGAAGGATGTGAACATGGTGGAGTTTGTCTGA
 GCCACAAGGGATCAGGCTCCATATAGGACTGAGGACATCCAGCTAGATTACAAAAACAACATCCTGAAGG
 AGAGGGCTGAGCTGGCCCATAGTCTCTGCCTGCCAAGGATGTGGATCTGGATAAAGAGTTCAGGAAGGA
 GTACTGCAATAAATGGACCCTGAGCTTCTGGCTGGGCCAGCAGCTCTGTATCAAAGGACATCTCCCTGA
 CCCTCCTGCGGTATTCTGGCTTCTCAGCGGCTGGTCCGACTTACCTAGAACTTGGCCTAAACTTGG
 CAGAGCAGCTGCCTGTACTTTGCCCTTCTAGAAATCGCCACCCCTCATCTTGGTGAGCAACTGTGGGTTT
 CCTAGAGACTCTGGTATAGTACGATTGCTGCCCTTCAGTCACCTGTGCCCACTGATGGTCTGACCCCAA
 CTAAACACAACAAGATCCCTTGTTAATATCCACCAATGCAAAGTCCCTCGTGGCCTCTTACTGCTAG
 GGTGAGGAAGACACTTAAAAATCCAGGTAAGACTCCCTAGCCACCAGTTAAACACATTAGCCATTGTCC
 TGGGGGGGGTCTTCTGAGCTGCATCGTGCCTGTGTACTGCTCAGAGCCCTGCTGTTATAGGTTCTGAC
 TCATGGGGCCGCTTGTGCTTTGGGCACTTGAGGCTAGCCAGGGCCCTTTCTGCTTCTGATTCTCT
 TTCTGCCCAATGCCTCCCAAGAGCTACACCAGCAGTTACTGGGTACCGTATGACCCTTGGCCTTGACATC
 CCTCCCTAGGCTGGAGTCTGGGTTGGGGCCCCATTTGCTCTGTTTTGGCTGAAGATGGGGTGAAGAT
 TTGGCTGAGTGGCCTATGCTGTACATCAAACAGCTATCATTTACTCCTACTTGGGAAGTTTTCATGTGA
 CAATAAAGATACATCTGAAAAAAAAAAAAAAAAAAAA

Restriction Sites: RsrII-NotI
ACCN: NM_024263
Insert Size: 1329 bp

OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC026438 , AAH26438
RefSeq Size:	2207 bp
RefSeq ORF:	1329 bp
Locus ID:	74761
UniProt ID:	Q9DBV4
Cytogenetics:	4 87.58 cM

Gene Summary:

Transmembrane protein which can modulate activity of various signaling pathways, probably via binding to integrin ITGAV:ITGB3 (PubMed:18366072, PubMed:22492581, PubMed:29702220). Mediates heterophilic cell-cell interactions in vitro (PubMed:18366072). Inhibits osteoclastogenesis downstream of TNFSF11/RANKL and CSF1, where it may function by attenuating signaling via integrin ITGB3 and MAP kinase p38 (PubMed:22492581). Plays a role in cartilage formation where it promotes proliferation and maturation of growth plate chondrocytes (PubMed:29702220). Stimulates formation of primary cilia in chondrocytes (PubMed:29702220). Enhances expression of genes involved in the hedgehog signaling pathway in chondrocytes, including the hedgehog signaling molecule IHH; may also promote signaling via the PTHLH/PTHrP pathway (PubMed:29702220). Plays a role in angiogenesis where it suppresses migration of endothelial cells and also promotes their apoptosis (By similarity). Inhibits VEGF-induced activation of AKT and p38 MAP kinase in endothelial cells (By similarity). Also inhibits VTN (vitronectin)-mediated integrin ITGAV:ITGB3 signaling and activation of PTK2/FAK (By similarity). May play a role in the maturation and maintenance of the blood-brain barrier (PubMed:14603461).[UniProtKB/Swiss-Prot Function]