

Product datasheet for **MC203007**

LOC433886 (NM_001007590) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: LOC433886 (NM_001007590) Mouse Untagged Clone
Tag: Tag Free
Symbol: LOC433886
Synonyms: MGC
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC025881
AAGCTAGAGGTAGGAGAGGCTGTTTGCTGACTGAGTACAGATCTCAGTGTCTACCTAGGTCTGGAGGG
ATGGGTAAACAACTGCCAGGAGCAGGGGACTTTTTAGTGGTACTTTGATTTTTGCATTATTGCTTTT
TAGATGATAGAAATACTAGCCATCTTTATCCAGCCAAAGCTGCGGATAGGATGGTGTGGTTATCTGAGGT
TTGGGAAAAGCTTAATATTGCCCCAGTTCCTTACATCACTGAAGAAGCCGAAAGCATTCACTCCATCAAT
AGAATTTTGCATCTAGGCCTCCAGGGGCCAGAATTGATCTCAGAACTCCAGCTAGAAAGAAGCAGAGG
AGAAATGCATTCCATTATCGTAAACCCTGTCTATGCCCTGGGACTCTCCGGGACAAGGACAAACCTGA
CGGGAGCCCCAATGTCCTGTTTATCTCTATTTTCAGTCTGTCCCCTGCCCTGAAGGCCAACCCCTTAGGG
AAGGAGGCCCTGAGGCTGCAAGTTTACACAGGCTAGTAGTCCCGAAGGTGCCTTTTTTCAATGGATGGG
GGCCACAGCTAGTTTTAGGTGCCTGTGGTCAGTATTGCTGCATAGGACACCAGAATGGCTAAGGGAAC
CCAGCCAGATGTGACTCAGAAAAATAAACACTCTCACCTCCCAGAAGCTCTGGTTTGGACACACATTCA
ACAGGTACACTTCTTTGTGGCCATTCATATTTTCAGGTAAAGGGGGTTTATAGTAAACTTAGCACAGT
CTAGTAGGGCAAGAACTGCAGTGAGCCCCAGGGTCTGTGCTTTAGGGGAGACAGGATAATAGGGGTCTGT
AGCTTGAGCTCAATGTATACATCCAGATGGAGAGAGCCATTAAGCTGATGTTGTTATATATATGTAAC
TGACAACAATATCCCCACTTCTGCTGTCTACCCAACATGGTCTTCTACTCTGTAGACATTAAGAGA
GTCCAGTCACTAAAAAAAAAAAAAAAAAAAA

Restriction Sites: RsrII-NotI

ACCN: NM_001007590

Insert Size: 336 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC025881, AAH25881</u>
RefSeq Size:	1010 bp
RefSeq ORF:	336 bp
Locus ID:	433886