

Product datasheet for MC203012

2310011J03Rik (NM_025521) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 2310011J03Rik (NM_025521) Mouse Untagged Clone
Tag: Tag Free
Symbol: 2310011J03Rik
Synonyms: AI452186; AW545358
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >BC023429

```
GTTCTTGGAAACGAAGATAGGCAGCGGGGTGTTCCAAGAGCAAACCTGGCGCGGAGATCGAGCCAACTCCG
AGGTTCCGGGAGGAGGGAAAGCGAGACGGGGTAGGGACCTCGGCGCAGGTCTGTAAGGGGTGCGTAGGAG
TTCGCGAGGTGGCGGGCCGAGCCGAAAGTGGACCCCAAGTTAGCGCGCGCTCCACCTACAAGATGAACT
CCAAGGCGAAGAAACGCGTGGTACTGCCACTCGCCCCGCGCCTCCACGGTGGAGCAGATCCTGGAGGA
CGTGCGAGGGGACCCGCTCAAGATCCTGTTTTCACTGCCTGGCCCCGGAAGACCCCCAGACCCATCC
CCAAGGGCTGAGGACTCTGAGATCCAGCAGGAGCAGATCTACCAACAGAGCCGGGCTACGTGGCTATGA
ATGAGAGGCTGCGGCAGGCGGGGAAGCTCTGAGACAGAAGTTCGATGGCCTTCGGCAAGCTGGTCAGCG
GCTGGAACAGGATATCAGCCAGGTGACTTACGCCACCTCTAGAGCCGAAGTTACACCTGATGTCCTTGG
TCCAGGCTGGGCACCTAGCAACCAGAGGCTCTGGGTCAACATTTGCTGTCCCAACCCTGGCCCCGGGGAG
ACCCAGTGGCCACCTTGCCAGCTAGCAGCTTTGTTCCCTCTTCCATCCTTAATGGCCTTAATTCATGA
TGGGACTTAGGGACATGGCCACCTCCTCATTATGGGGACACTGGCCTAAATTGCATCAAGGACTTGATTG
CAGGTCCGGGCTCAAAGCCAGTCTGCCCATAGTGACCGGTGTGCTCTTCTGTCCATGTACCCCTGGGC
AGTGAAGGTCTTCATCCCCAGCATCCACAGCAGGCCAGCCACTGGCATCTTGGTTGAGAGAATGGGCTAG
CTGCACGCTTCATAGAGCCTGGCTCAGGCCCTCTAGGGCTACACTGCCTTCTCAGTTCTCCCTGGGGCCG
GGGAGACTCATCAAGGGAGTTTCCATTTGGGGAGCCAGGGTCCAAGAGGAACCAAAGGCCAGGTCC
CTGTGTGGTGAATGTTGCCCTTCTGGATCTGGGGTTAGTTTGCCTTCTGCAGAAGACCCACCACT
CTTGGGGTGGCTGGACTTTAGTTCTGGAAGATTAAGATGTTTGTGATGCCTGAACACAGTTCTCCCTTAAG
CCAGGCTAGATCACTGGGTACCCAAGCCCCCTGCCACCCCTTGGCAGCTCAGGCGTAGTACCACACT
GCCTCCTTGGGTCCAAGCTGGCTGGAAGTGTGTTTGTGAGAAGTGGCCTTAGCTTGTGGCTACATCATT
CTGACACTCCACATTCCTGGCGTCAAGGGTCCAACCTTGCCCTGACCAGTGCAGCTTGGACTCCTGGTC
CCCACCTCCACATGACACCTCCACCACAAATGACCAGTGGCTGATGGCCAGGTACAGCTTGGCTGAACC
TCAGTCAACCTCGGCTCCAGCATGTGGGATGTGACACGCTGTCGAGTCAAGTGGGAGCTCAGAAACACC
TTCTGGGGTTGAACAGAGCTGTAGAAGAAGTCTTGAACGACCCCCCCCCCAGAGCCAGTGATAAAG
ATGAAGCTACAGATCGAAAAAAAAAAAAAAAAA
```

Restriction Sites: RsrII-NotI



[View online >](#)

ACCN:	NM_025521
Insert Size:	330 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC023429 , AAH23429
RefSeq Size:	1641 bp
RefSeq ORF:	330 bp
Locus ID:	66374
UniProt ID:	Q9D7E4
Cytogenetics:	10 C1