

Product datasheet for MC203021

Cdpfl (NM_197998) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Cdpfl
Synonyms:	2210021J22Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC025858
 GGGGGTCCCGAGGGGAGGTGACAGCTGCTGTTGCCATCAGACCACTGACCCAGCAATTTACGGAGAGGC
 GGGGCCTCCAAGGCCTTGAATATAGAAGTGCCATCAAACAGTTAAGAATGCAGCCTTATCAAGTGACTA
 GGATACACCTCCCAACCTCCTGATGGCATCTGAGACAGAGCCCGCCCTGGGCACGTTTGAG
 TGCCAGCTTTGTCCTTATCAGCTCCTTACAGCTATGTGGGCGAGAAGCCCCAGACACCCAGGCTGTTG
 TCCTCCTGGAGGAGAGCTACATCATGAAAGACCCCTTCTCCTCAGACAAGGCCAGGTTTCTGGTCTCGG
 CTCTCGGTGCACTGTGTGACGAGGCTGGTGTGTGTGGGCCCGGACTGCAGCCTGTTCTACTCCAAGAGG
 GTCTGCCTTCCCTGTGTCCAGGAGAATGAGTGCTTTCCCTCAGGAAATCCAGCAGGATGTAGAGAAGA
 GAAAGAGCGCCAGCAAGAAGCACTCCAGCCGTCCTGACTGGCAGGGGCACAGCCCTCCTGCGAAGCATT
 GGGCTGGGCTCAGCCAGGATTGAGGCTGGGAGCTGTGGACCGGTCTTTGAACTGGGTGGATACCTC
 CAGCCTCGGAAGCTGAGGAAGTGGCAGGGCTGAGGACAGACATCAGGGGTGTATCCAGCTGGAGCTCTA
 GGTCTGGCCTCTGGAGGCCAGGAGAGATGGAGAGGGACCTCCTTTGATTGACAGCTGGTAGACACTCCCT
 GTAGATGCACAGTGGGAAGACACCTAGCTTCTCTCCAGACTTCAAACAACACAATGTCATCTGTCTAG
 GGAATTCAGCCTCCGAGGGCCAGAGGTAAGCAAGGGTGACGTACATCCCTGTAACCTCTGGAGATGA
 TTAAACCTTGACTCTGTATCATATTTATCATTAAAAACCTGAAGGATAAAAAAAAAAAAAA

Restriction Sites: RsrII-NotI

ACCN: NM_197998

Insert Size: 360 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC025858</u> , <u>AAH25858</u>
RefSeq Size:	973 bp
RefSeq ORF:	360 bp
Locus ID:	72355
UniProt ID:	<u>Q8R3A2</u>
Cytogenetics:	15 E2