

Product datasheet for MC203027

Flad1 (NM_177041) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Flad1 (NM_177041) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Flad1
Synonyms:	A930017E24Rik; Pp591
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC025817

```
TGGGATGATACGAGATCAGTGATCTGAAGCCTGGTTCCTCAACCCCGAACTCAGGGTATGGCCTCCTG
CATCCTGTAGACCTGGCGCCTCTCAGTACATACGGCCCAATTTGGGGTCTGAGTGGCTACTGGAAG
CCTTGCAAAGAGCACAGAGGGCAGGACCATGGCATCTAGGGCCTCTGAAGTGCCTCCAGGATCAGGGCG
CAGCGTGACAGCAGGCATCATCATTGTTGGAGATGAGATCCTTAAGGGACACACGCAAGATACCAACT
TACTTTCTTTGCCGACCTGCGCTCCCTGGGGTCCAGGTTTGGCGAGTCTCTGTTGTGCCTGATGAAG
TAGCCACCATTGCATCCGAGGTCAACTATTCTCCAGGCGCTTCACTCACGTCTCACCGGGGGGCAT
CGGCCCACTCATGATGATGTTACCTTTGAGGCAGTGGCAGGCTTTTGGGGAGGAGCTCAAACACAC
CCTGAGCTCCAAGCAGCCATCAAACCTAGGAGGGGAGGTTGGGAGAAGCTGTCGATGGTGCCTCCT
CTGCCCGCTGCATTACGGTACAGATCCTCGAACCGCCACCCCTCAGATTCCCCTGGTCTCTGTCCG
AAATGTCTACCTCTTCCAGGAATTCGGAGCTTCTGCGACGAGTTTTGGAGGGACTGAAAGGGCTGTT
CAAACACAGCTGTGCAGTTCACCTGAAGGAGTTGTATGTGGCCGCTCAGAGGGCTATCGCCCCCA
TCCTGTGGAGGCTCAGGCCACTTTGGGCGTAGGCTTAGCCTGGGTTCTACCTGACTGGAGCAGCAA
CTACTTTCAAGTGAAGCTGATTCTAGACTCCGAGGAAAAAGGCCCTGGAGGAATGCTGGCCTATCTG
ACTGCAGCTGTGCCCAAGGATCACTGGTCCCCTACCAGCCCGATGCTGTTGAGAAGGCTGGTGAGGCTG
TGTACAAACTCGCGGAGTCAGGCTTGTCTGGGAAGAAGGTAGCAGGTGCCCTACAGACCATCGAGAC
TGCCCTGGCTCAGTACCATCTCTCCAGCTCTGCGTGGGCTTCAATGGGGCAAGGACTGCACCGCACTC
CTGCACCTTCCATGCAGCCGTCAGAGAAAAATTCCTGATGTTCCAAAACCTCAGATCCTGTACA
TTCGTAGCATCTCGCCCTCCCTGAGTTGGAACAGTTTTCTGAGGACACCATCAAGAGGTAATCTCCA
GGTGTTAGAAGCCGAGGGTAACATGAAGCAGGCCCTGGGTGAGCTTCAAGAAAAGCACCCCAAGCTGGAG
GCTGTCTGATGGGCACTCGACGACTGACCCCTACTCCTGTAGCCTCAGCCATTTAGTCCCACCGACC
CTGGCTGGCCTTCAATCATGAGAATCAACCCGCTGCTGGACTGGACCTACAGAAACATCTGGGAGTTTCT
GCGGCAGCTGTTGTCCATACTGCATCCTATATGACAGAGGGTACACATCCCTGGGGAGTCGGGAAAAAC
ACAACGCAGAATCCGGCCCTGAAGTGCCTGAGCCAGGAGGACACCCTGTCTACCGTCCAGCATACTAC
TGGAGAATGAGGATGAGGAGAGAAAATCCAGAATGTGATCTCACACAAGCTGGGGACAAGACAGAGGGAA
GGGGCTGGCCTGGGCTACCTGCATGTCAGTAAAGTGTACCTCACACCCCAAAAAAAAAAAAAA
```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_177041
Insert Size:	1479 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC025817 , AAH25817
RefSeq Size:	1746 bp
RefSeq ORF:	1479 bp
Locus ID:	319945
UniProt ID:	Q8R123
Cytogenetics:	3 F1
Gene Summary:	Catalyzes the adenylation of flavin mononucleotide (FMN) to form flavin adenine dinucleotide (FAD) coenzyme.[UniProtKB/Swiss-Prot Function]