

Product datasheet for MC203036

Rmdn2 (NM_201361) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rmdn2 (NM_201361) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Rmdn2
Synonyms:	AW061290; Fam82a1; mRMD-2; RMD-2
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC024059

```

CCACGCGTCCGCCGTTATCTGGGTGACACACGAAGGCCACCACGCAGGACCTGGGACCCGACGGCTGCA
CAGGGTCCACTCGGCCATCCTCTCCTGGCTCGCCGCTCGCCCGGGTCCCGGTTCCCTGGCTCCAGC
CCCAGAGAAGCCTGCGCAACGACTCCGGAAGTAAAGGCCAGAATGCCCAATCCACAAACAAGAGCTA
ATCCTTGGCATCATGGCGGGGACTGCTGGGATCAGCTTGTGGCCTTCTGGTACCACAAGTCTCAAGC
CAAGGACGACCATGAACTTCTAAGCTCTTTCTCTGGGTAAGAAAATTTGGCTCACTAACTCTGCCAGA
GGAAAGCCACAGTGCACAAGGAGCATCGGTGGTCTTTCAAAGGAGGCAACTCCAGATCCTGGAAAAGCTA
AACGAGCTACTGACCAACATGGAAGAACTGAAGGAGGAAATCAGATTCTTAAAGAGACAATCCAAAGC
TGGAGGAATGTACAAGATGAATTCGGAGGGAAGGTGACAGTTCATAAGATCAGTCTCAGCACAGAGC
GAGAAAAAGAAGGGGACAACCGTTCAACGCTCAGCAACGAGTAAACGCTCAGAGGAGGCAGAAAAGTGA
GGGGGTATATTACAGCTAACACGGACACAGAAGAACAAGAACTTCCCGTCCCAAGGCTCTTAACACAC
ACATAGAGGAGCTAAAGTTGGATGTCCTCTCAGAAAGGACAGCATTTCGCGACGAATGAGTCTCACAA
AATGGAGAGTTTTGAGCTACTCTGTGACCACAAGGAGAAGTTTAGTGAAGAAACTGAGTTTCTGTGGCGC
CTGGCCCGAGCCTATGGAGACATGTATGACCTATCCACAAGTACACAAGAAAAGAAGCATTATGCTAATA
TTGGAAAACCTTAGGAGAGAGAGCTATTACCAGAGCGCCATGAATGGACACTGCCATCTGTGGTATGC
AGTTTTGTGTGGTTATGTATCTGAGTTGAGGGCTTGCAAAACAAAATTAAGTGTGGACATCTATTCAAG
CATCTAGAAATAGCCATCCAGCTTTTGCCTGAAGAACCCTTATTATATTACCTCAAGGGAAGATACTGTT
AATCTGTCTCAAGACTCAGCTGGATTGAGAAAAAATGGCTGCTACTCTGTTCCGAGAAAATACCATATTC
AACGGTTCATGAAGCTTTGCACAACCTTCTTAAGACTGAAGAAGTGAACCTGGTTATTCCATGTCCAAC
TATATGTATACAGCCAAGTGTATGTTGAACTTGGGAACTCAGGAGGCTTGTAAATTCGTAACTTGG
CATTGCTACTGCCTGTTGTACCAAAAGAGGACAAGGATGCCCATAAAGAGGTGAAGAAAATGATCAGTTC
TTTGAAGAGGTAATAAAAGACCTTATTCTCAGCAAACCAATACCGTCTTCTAGAATTTAACAGATTG
TACTTTTATAATGCCTTCTCCTTAATGTAGAGATTATGACTCTCATTGATTGAAAGGCTAGCCTTG
CCTTCTAGACCCACTGCACTATGAGACAAATCATGTTCTTTGGAGGCAGTACTTACTCTGATGTTAGTAT
GTAACGATCTCTGGACATGTATACTTCTTTGTTTATTTCATGTCAATATACTACAGACTTTAGAATGTTTG
TGTAACCAATACGTAATCTAATAAAAATGCACCCTAGCAAATGTTATTGTTTCACTAGAAAAAATAAAAAA AAAA

```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_201361
Insert Size:	1191 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC024059 , AAH24059
RefSeq Size:	1754 bp
RefSeq ORF:	1191 bp
Locus ID:	381110
UniProt ID:	Q8BSE0
Cytogenetics:	17 E3