

Product datasheet for MC203053

Arxes1 (NM_029541) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Arxes1 (NM_029541) Mouse Untagged Clone
Tag: Tag Free
Symbol: Arxes1
Synonyms: 6530401D17Rik; Spc22/23; Spcs3
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >BC030351

```
CCACGCGTCCGGCCAGAGGAAGGCAGAGGGACCGCGGAGCGAACTGCACGCCCGTGCAGAGCCAGCCAG
CCAGCCAGCCAGCCAGCGAGCGCGGTGGAGCCCGCTGGAGGGAGCCCGGATTGAACACCTGCAAGA
TGAAGTCCCTGTGTCCCGGGCCAAGTCTGCTGTTGCGCTTCACGCTGAGCGTCATGGCGGCGCTCACCTT
GGGCTGCATCCTCACCACCGCCTTCAAAGACAGGAGCGCGCCCGTGCCTGCACGCTCTCCAGGATCCTG
CTCAAAAAAGTGGAAAGACTTCACCGGACCCAGAAAAGAAGAGCGACCTGGGCTTCATCACCTTCCACATCT
CTGCGGATCTGGAGAAGACTTTCGACTGGAACGTCAAGCAGCTCTTCTTTATCTGTGCGCAGAATATTC
CACCAAAAGCAACGCTGTGAACAGGTGGTCTCTGGGACAAGATCCTTCTGCGAGGGCAGAAACCCCAAG
CTGAACCTGAAAGACGTCAAAGCAAGTATTTTTCTTTGACGATGGACACGGGCTCAAGGGAAACAGGA
ATGTCACCTTGTACTCTCTCTGGCAAGTGATACCGATTGCTGGGATCCTGCCTCTTGTGACTGGATCTGG
ACGCGTGTCTGTCCCGTTCCAGATTCCATGAAATAGCCACGACTTTTTGAAGAATTCTCTCTGCATTG
GAAACGGGGACATCTGTGGACTCAAGAAATCGCATCTCCTCATCATTTCCCTCAAGTCTTGTGACTTCT
CTTTTGTTTTGTATTTTAAAAACATGCTGTTTTTGTCTGGGACTGGCTGCACTCATCCTCATGCTTCC
CTGAGCCATAACATAGAGACTATATTCATTGATCTTAGTGTGTGTGTGACATTGGTTATCACTCTTCTCT
CTGTATAAGAATTAACATCAAAATTTATAGTTGGAGGGAAAGGCTAATAGAGTTCTGGAACCTTCTTTAT
TCTACCAATGAAACAACAAGAAGTTCACCAGAGTAGGTTATTATTAGACAGTTGGATAATGAGGATTGA
AAATAATACACAGTTTTTGCATATTCACCTACATTATCCCATGTGATATATAAAATTTGGGAGAATTTTA
AAATTAGATGTTTTGGTATTGTCAACAATAACTGATTTTTATTCTGAGATTCTTGTAGGTCCACATTGA
ATGCCGTTTTCTATATCTTAGGGTGTGTTGATTACTAGAAAAATTTGATTTTGTGTTTTCGTAAGAGTC
TGACATGCAAAATAGTTAAACACGAGTCAGAAATCTTACTGTTTCTTCTTCCAGCAAAATGTGACCA
TTTAAATAGTGAACGTATCTTGAGGCTGTAATGATTCTGCCTCTCTTTAAGTGTGTGCTTACAGTGGT
TTCTCAACAAAGATTGCTGTGCAAGTGTATCTTTGTCTTAGCACTGTTGCTGAAGTTTAAAAATGTTTGA
TGTTTAAACATGAGATCGTATGTTCTAGTATGTTTAAATGCCATTAATAAATGCTAACTGCTTAAAAAAA
AAAAAAAAAAAAAAAAAAAA
```

Restriction Sites: RsrII-NotI



[View online »](#)

ACCN:	NM_029541
Insert Size:	1245 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC030351</u>
RefSeq Size:	1558 bp
RefSeq ORF:	1245 bp
Locus ID:	76219
UniProt ID:	<u>C0HK79</u>
Cytogenetics:	X F1