

Product datasheet for **MC203154**

Gstm6 (NM_008184) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gstm6 (NM_008184) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Gstm6
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	<p>>BC031818</p> <p>CAGACTCATCAGCTGAGCATTACAGAGTCTTTCCCGACCAGTGACTGACTGCACAGTCTTCGAGGCTGGA ATAGCAATCATGCCCGTACTCTGGGTTATTGGGATATCCGTGGACTGGGTCACGCCATCCGCCTGCTCC TGGAAATACACAGAAACAGGCTATGAAGAGAAGAGATACGCCATGGGGGACGCTCCTGACTATGACAGAAG CCAAGTGGCTGAATGACAAGTTCAAGCTGGACCTGGACTTTCCCAATCTGCCCTACTTAATTGATGGGTCA CACAAGGTCACCCAGAGCAATGCCATCCTGCGCTACCTTGGCCGGAAGCACAACCTGTGTGGAGAGACAG AGGAGGAGAGGATCCGTGTGGACATTTTGGAGAAACAGGTCATGGACACTCGAATTCAGATGGGCATGCT TTGCTACAGCGCTGACTTTGAGAAACGGAAGCCAGAGTCTTGAAGGGCCTCCAGATCAGCTGAAACTC TACTCGGAGTTCTGGGGAAGCAGCCATGGTTTGCAGGGGACAAGATCACCTTTCAGACTTCTCTCGTCT ATGATGTCCTTGATCAGCATCGAATGTTTGAACCCACGTGCCTGGACGCCTTCCCAAACCTGAAGGACTT CATGGCCCGCTTTGAGGGCCTGAGGAAGATCTCTGCCTATATGAAGACCAGCCGCTTCTTCCAAGTCCT GTGTACTTAAAAACAGGCCACGTGGGGCAATGAGTAAGACCATGCATGGAGTGCCGGTGTGGTGGGGAAG ATCTATCAGTGTGCCTGTCACCGTGAACACTGACCGGGCCACAGCTAGCTTCCAGTAACCTCTCTTCATC AATCTTCCATTTTTTCCCTCTGTCTTTTCATTAACCTTCCCTCTACAAGAACTCCTTGTCTCCCTTTT CACTCAAGCCTCCAATGTCAGCTCCCTGCTCCTCAGCAAAGTCTTTGCTTCCCTTGTCTTCTTCTGCCAT GTTTCTAACCAGACAATCTTCACTGCACTTTGAGCAGCTTAGCCCAAAGATCCATACTCCCAAGTATGC AGGCTCCATCTCCAGAGTGGGAAGGCCAGTCTTGTGTGAGCTCTAGTAAAGTGTATCCATACTAAAAA AAAAAAAA</p>
Restriction Sites:	RsrII-NotI
ACCN:	NM_008184
Insert Size:	657 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC031818, AAH31818</u>
RefSeq Size:	1129 bp
RefSeq ORF:	657 bp
Locus ID:	14867
UniProt ID:	<u>O35660</u>
Cytogenetics:	3 F2.3
Gene Summary:	Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.[UniProtKB/Swiss-Prot Function]