

## Product datasheet for MC203165

### Pole4 (NM\_025882) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Pole4 (NM_025882) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Pole4
Synonyms:	2400007P05Rik; 5830430F06Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC023189

```

GGGGATGGCGGCAGCGGGGCTGCTGGGAGCGGGACGCCAGAGAGGAAGAGGCTCCCGGAGGGGAGGCA
GCGGCTTCGCAGGCCAAGCCCCGACGAGTGCACCCGGGGAGTTGCGCTCTCGAGGCTGCCTCTGGCGA
GAGTAAAGGCCTTGGTGAAGGCAGACCCTGACGTAACGCTGGCGGGACAGGAAGCCATCTTTATCCTGGC
GCGTGCCGCGGAAGTGTTCGTGGAGACTATCGCAAAAGATGCCTACTGCTGTGCTCAACAAGGAAAGAGG
AAAAGTCTCCAGAGGAGAGATTTGGATAATGCAATAGAAGCTGTGGATGAATTCGCTTTTCTAGAAGGCA
CTTTGGATTAATTGTGGAGCTGGATGGTCTGAGCCTTGGCTGCATCTGTACCTTTCCCTCTGTAAGCC
TCAGCTTTGGAGAAAGTCTGCGCATGTACACAATCTTCTACTCTTTTTACCCAAGATAGAAGAAGCAGA
GAGCTCTGGATAACGTTCTCCACCAAGTCTCCATTACTTGAGTCTGTCTCAATTCAAGCCTGGATG
TAGCTGATGCTGTTTGGGCAGAGATGAAAGACTGTTCTGCATAGCTGGCTTAATGAGCACAGACTAGA
ATAAAAGTCTTTGTTTCAGCTTCTTTCCTGTGTTTTATATTGCCTTTGCTCTTGCCTGTGACTGGGAAC
CCATGTGATCAGGCACAGACTTTCTCATCTTTCTACCAACCACTCCAAGTCAGTCTTTTTATTAATTT
GTTCTTTCTCAGCATGTATTGGAATGTGGATCGATAATATATGGGATACCTCTTCTCCTTAGAAGTCAA
ATCGATGTGGTCCCAACAGTTTTATTCAAGTATAATTCATGTGTCATATTCATTCAAAGTAAACATTC
GGTGGTTTTTAAAAAATATTCATAGATATATGTAGTGAGATTTTGAACATTTTCTGTGAATGGATTCA
CATGAGTCTTGTGTATGGCAAAATTTAGATAATGTTATTTAAAGTTAATCTATGCTATGATATGATTTA
TACCCTTAAAGTTGTGAAATATATGTGATATAAAATTTCCCATTTTAACTGCTTTTCAAGTGTGTTGCTCA
GTGAGATTGAATATTTTCAGAGTGTGACACATGGAAACAGTATCTTCTTTTCTCATCTCAGTTAGGCTAA
TTCTCAGCCTGGAACATTGTGCACATGGAAACAGTATCTTCTTTTCTCATCTCAGTTAGGCTAA
ATTAATTTTGTGTTGTTGTTATTTTGTGAAGTGTGGGCAGTGCACCTTTAGTACCTAGAATGTAGTGCT
TCAGGAAAGCAGCAGAAAGCTGTGGCTAGGGCTTCTACAGTGTCTGCAGTGTAGATCTCTTAGGATAG
AGTCTGCAGGATCTTCTCCCTAGAAATGATACATTCAGTCTGTAATACGTTAAAAATAAACTTTTCAGAAA
AACTACTACATTTGTAGAAAAGTGTGTTTGGTCTTACTAAGTCCACCCAAAAGTACATTTTCAGGG
TTACCAGAAAAGTACTGACATTTATGGCTTATTGAGTTTCAAGTTTTTATTTGAGTTTCAAGTATTTTTTCCA
TTAATATCCTTTGTTTGTGAGATCTAATCGAGATTATAAATTGTTTTGCTTCTTGAAGAAAAA A

```

Restriction Sites: RsrII-NotI



[View online »](#)

<b>ACCN:</b>	NM_025882
<b>Insert Size:</b>	357 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC023189</a> , <a href="#">AAH23189</a>
<b>RefSeq Size:</b>	1681 bp
<b>RefSeq ORF:</b>	357 bp
<b>Locus ID:</b>	66979
<b>UniProt ID:</b>	<a href="#">Q9CQ36</a>
<b>Cytogenetics:</b>	6 C3
<b>Gene Summary:</b>	May play a role in allowing polymerase epsilon to carry out its replication and/or repair function.[UniProtKB/Swiss-Prot Function]