

Product datasheet for **MC203195**

Mrps11 (NM_026498) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Mrps11
Synonyms:	0710005I03Rik; C79873; MRP-S11; S11mt
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC028279
 CCACGCGTCCGATTTCAGTCATGCAGGTTTAAAGAACTCCGGGTCGTGGCTCCTGCTGCTCGGCTGG
 ACATGGAATGACCAAGTTTGTGGCATTGCCCCAGCTAAGAATATCCACACCGGAGCCCCGAGATTAGAA
 GACTCAGCTGCCAGGCAGAACTGAAAGGGAAGCGGCTCCAAGTCGCTTCAGCCTTTACCTCCGGTTC
 CAGGACAGGAGAGCTCTCTGCAGTGGGCAGGGATGAAGTTTGAGGACGTCCCAATTGCACACATTAAGC
 AACCTACAACAACACAGATCCAGGTAGTCTCTGCCACCAATGCGTCCCTGGCCCGCGCTCCTGTGGC
 ACAGAAGGCTTTGCAATGCCAAGAAGGGCACCAGCATCGCCGCACAGACAGCAGGCATAGCTGCCGAG
 CAAAAGCCACAGGAAAGGGAGTGACCCACATCAGAGTCGTAGTGAAAGGCATGGGGCCAGGACGCTGGT
 TGCCATCAAAGGACTGACCATGGGGGTCTGGAAGTGATCTCAATCACAGACAACACCCCGTCCACAC
 AATGGCTGCCGACCCGGAAGGCTCGAAGGCTGTGACGGGAAGGGAGCCTGCACTTGACCTGGCCTCAG
 GTCTTGACTCCTGTGGATCCTACAGAAGGCCACTGTGGAAGCCGCTGGAGCCGCTGTTCTCGAGCCCT
 CCTTGGGGAACACTCCTCTCCTCGCACCAGGCTGGAAGTGCCAAGCTGGCCTCTTGACATGAGGGAG
 TGGGAACCCAGGCCACTGTGGCCAGCTCAGCCTCTGCTCCAAAATAAAAACTATTCTCTGTAAAAA AAAAAAAAAA

Restriction Sites: RsrII-NotI

ACCN: NM_026498

Insert Size: 576 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC028279</u> , <u>AAH28279</u>
RefSeq Size:	852 bp
RefSeq ORF:	576 bp
Locus ID:	67994
Cytogenetics:	7 D2