

Product datasheet for MC203348

Olfr722 (NM_146494) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Olfr722 (NM_146494) Mouse Untagged Clone
Tag: Tag Free
Symbol: Olfr722
Synonyms: MOR241-3
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >BC051403

```
CTACTTTTTTAAAAATATTGGACTGCTACTAAGGAAGTGATGGAGATCAAGAACAGCTCAGTAGTAACAG
AATTCATCCTACTTGGCCTGACCCAGTCTCAAGAAGCTCAACTCCTGGTCTTCGCACTGATTTTCAGTTTT
CTACCTAATAATTCTCCCTGAAAATTTTCTATCATTTTTCAACCATTAGATCAGACTCTGGACTCACAGCC
CCACTTTACTTCTTTCTGGGAACTTGGCATTCTGGATGCCTCCTACTCTTTTCATTGTAGCTCCTAGAA
TGCTTGTGGATTCTTCTGTGAGAAAAAGTCATTTCTTATAAAGCCTGCATCACTCAGCTATTCTTTCT
GCCTTTCTTGGAGCAGGAGAAATGTTTCTTGTGTGATGGCCGTTGATCGTTATATCGCTATATGC
CGTCTTTATACTACTCCACTCTCATGAACCTAGAGTCTGTTATGCCTTATTGTTGGCTCTGTGGCTTG
GGGTTTTGTCTATTCCATTGTACAAGTGGCTCTTATCCTAAACTTGCCTTTTGTGGGCCAAACCAGCT
GGATAACTTTTTCTGTGATGTTCCACAGGTTGAAAGCTAGCCTGCACTGATACCTTTGCTGTGGAGCTC
CTGATGGTTTCTAACAGTGGTCTACTTACCCTCTTATGCTTCTTGGGACTTCTGGCCTCCTATGCTGTCA
TTCTCTACCACGTAAAAGGACATTCTTCTGAAGGGAAGAGCAAAGCCATCTCTACATGCACCACCCATAT
TATAATTGTGTTTCTCATGTTTGGGCTGCCATTTTTATCTATACCCGCCCTTCCAGGCTCTCCAAGCT
GACAAAGTTGTTTCCCTTTTCCACACGGTCATCTTTTCTTGTGATGAATCCTGTGATCTATACACTTCGCA
ACCAGGAAGTAAAACTTCTATGAGGAAGTTGTTAAGTCAGTATGTGATTTGCTGAGTGGATAAAAGAAA
GAGAAAGCAAAGAGAAAGAGAGTGAATAACAGTTTGAATTAATAAATGCAATTATTCATGCATTG
AATCAGTTCTTTAGCAAATGCTGTGTTAGTACTTCCAGTTTTCCGACTTTATTTAAGAACACAAAATGA
GACAGGGTTGTATTCCAGTGAATGAAACCAAGCCAAACCTTCAACTCGAAAAATCACTAGAATCAATG
CTCTGGTATTTTTAGACCTTAGCTATTGAGTAAGACACTTGAGTAGACTTTAACTGAAATCAAGAGTAT
TAGCAGAAAACAAAATAACAGTTATGAAACAAGCATCATAGAAAGCACTGGAACCTAAAAGATTTGAAGCTTT
GAGTACAATTGAAGGACATAAAAAGGAAGAGAATAGTTACAATAGAGTTAAGGGATGTAGAAGCAGCAGAA
TGTGACCTGGATAAAGATTTTAGAGCATTTTTTAGGTAGCTTTAGTTATTGAAATCCTAAAAATCATACAA
CTGAAATTTATAATGATTATCCAGTGTCTTCTATTTGATTTAAGACAAATACTAAAAATCTGTATAAAA
GAAAAATATATGATTAACAAAATAAAATTAATGGAATGTAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
```

Restriction Sites: RsrII-NotI



[View online »](#)

ACCN:	NM_146494
Insert Size:	927 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC051403 , AAH51403
RefSeq Size:	1655 bp
RefSeq ORF:	927 bp
Locus ID:	258487
Cytogenetics:	14 C1
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]