

Product datasheet for **MC203388**

Fam214b (NM_172691) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fam214b (NM_172691) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Fam214b
Synonyms:	B230312A22Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC071191</u> , <u>AAH71191</u>
RefSeq Size:	2811 bp
RefSeq ORF:	1617 bp
Locus ID:	230088
UniProt ID:	<u>Q8BR27</u>
Cytogenetics:	4 A5