

## Product datasheet for **MC203503**

### Pmm2 (NM\_016881) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Pmm2 (NM_016881) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Pmm2
Synonyms:	A1585868; C86848
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC046325

```

GTTTCTGACTGGAGAAATGGCCACTCTCTGTCTCTTCGACATGGATGGGACCCTGACTGCCCCGCGGCAG
AAAATTACAGAAGAAATGGATGGCTTTCTACAAAAATTGAGGCAGAAGACCAAAATGGAGTGGTAGGTG
GGTCAGATTTTGAGAACTGCAAGAGCAACTGGGAAATGATGTGGTTGAGAAATACGATTATGTGTTTCC
AGAGAATGGCTTGGTAGCGTACAAAGATGGGAAGCTATTGTGTAAGCAGAATATTCAAGGCATCTGGGG
GAGGATGTGATCCAAGACCTGATCAACTACTGTCTGAGCTACATTGCAAACATTAATACTCCCTAAGAAAA
GGGGAACTTTCATTGAATCCGAAATGGCATGTTGAATGTGTCCCAATTGGAAGAAGCTGCAGCCAAGA
AGAACGAATTGAATCTACGAACCGATAAAAAAGAACATATACGACAAAAGTTCGTAGCAGACCTGCGG
AAGGAGTTTGACGGAAAGGCCTCACGTTCTCCATAGGTGGCCAAATCAGCATTGACGTCTTTCCTGAAG
GATGGGATAAGCGGTATTGCCTGAGACACCTGGAACATGCTGGCTATAAGACCATTTATTTCTTTGGAGA
CAAGACTATGCCGGTGGCAATGACCATGAGATCTTACAGACCCAGAACTGTGGGCTACACAGTGACG
GCCCCGAAAGACACACGCAGGATCTGTGAGGGGCTCTTCCCTTGACCCGAGGAGTGCCGGACTGACCAAC
AACACTGAAAGGGGCATCTGATGATGAAAGAGAGCAGCCACCCAGAGTTACCCACCTGGGATAGAAGTGA
AGACTTCTGTGACATGCTCCGTCTGCAACCTGTGGTGGAGGCTGAGGCAGGAAGGAACACCACTATTCTG
TCAAAGGGTAAATGTTTCTGTTTCAAGCCAGCTTCCCTCTGTGAGCCACAGACGACAGTGTGGCCTGC
TGCCTGACAGATAGTGACCAAGTGACAACCTGATCCACAGTTTACACTTGCCTGGCTGGCTGGCCAT
CTCTATTTGCTATATCTTAAAGCCAGACTGCTTTCCTGTACCAGGAGGGGAGGACAGAGGACCCAGC
AAAGCAGCACCTCACTGCAGTCTCTTTAACCCAGGCTGACCAACCTAGGTAGCAGACGCAGTACCTA
AGGCCACAGCAGCAGGGGTGCTCTCCAGCCTCTTCTGGAGAGTGACGTGTCTCCAGAAAGCTCGGG
GACAGCCTTCCAAGAACACTTCCAGGCAGCTTCCAGTGCCCAACCACATATGTCTCCACTGCCCTAC
CCTTGGCTGCCCTACTCTGTCCACACCACAGCAGAATGCCAGGTCTCACCTGTGGCATCCTTGGGT
TTTATGGCAGTGTGTGAAGCCACTGAAGGGACAGCACCATACCTCACTTCTCAGTCTGGTCACTTTAA
CCCTCAGTGTCAAAGGGTGTCTACCCAGGTTGGCAGAGAAGTGGCTTCTGAGAGTGTCTTCCCCAC
TTGGCTCCAGTAATGCTGGAGCCCTAGCTGTCTGGCTGTAGCAGTTCCAAGCCAGCCAGATGCTGTG
TTCTACCTGCAGGTCCAGTGTCTCCACTGTGCCAGGTGGACTCCAGGCCTTCCAGTTCCTACAGTGGC
TATCACAGAACCCCTGCCCACTGTCTTCATTCAACTCTCCACCTAATAAAGGGCCTCTGCTCTCCA GGCA

```



[View online »](#)

<b>Restriction Sites:</b>	RsrII-NotI
<b>ACCN:</b>	NM_016881
<b>Insert Size:</b>	729 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC046325</a> , <a href="#">AAH46325</a>
<b>RefSeq Size:</b>	1754 bp
<b>RefSeq ORF:</b>	729 bp
<b>Locus ID:</b>	54128
<b>UniProt ID:</b>	<a href="#">Q9Z2M7</a>
<b>Cytogenetics:</b>	16 A1
<b>Gene Summary:</b>	Involved in the synthesis of the GDP-mannose and dolichol-phosphate-mannose required for a number of critical mannosyl transfer reactions.[UniProtKB/Swiss-Prot Function]