

Product datasheet for MC203515

Copz1 (NM_019817) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Copz1 (NM_019817) Mouse Untagged Clone
Tag: Tag Free
Symbol: Copz1
Synonyms: 5930435A22Rik; AA407760; D4ErtD360e
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC058524
GCAAGATGGAGGCGCTGATTTTGGAGCCCTCCCTGTACACTGTCAAGGCCATCCTGATTCTGGACAATGA
TGGAGACCGACTCTTCGCAAGTACTATGACGACACCTACCCAGTGTCAAGGAGCAAAGGCCCTTTGAG
AAGAACATTTTCAACAAGACCCATCGGACGGATAGTGAAATCGCTCTGTTGGAAGGATTGACAGTGGTCT
ATAAAAAGTAGCATCGATCTCTATTTCTATGTGATTGGCAGCTCCTATGAAAATGAGCTGATGCTTATGGC
TGTCTGAACTGCCTCTTCGATTCCCTGAGCCAGATGCTGAGGAAAAATGTAGAAAAGCGAGCTTTGCTG
GAGAACATGGAGGGCCTCTTCTGGCTGTGGATGAAATTGTAGATGGAGGGGTGATCCTAGAGAGCGACC
CCCAGCAAGTGGTGCACCGGGTGGCTTTGAGGGGTGAAGACGTCCCCCTTACAGAGCAGACCGTGTCTCA
GGTGTGCAGTCAGCCAAAGAACAGATCAAGTGGTCCCTCCTTCGGTCAAAGTGTCTGAGCCCAGTTCC
TTGCCCTCAAACCCACATCTGCTGTGACTTCTATCCAGGCCCACTTATGTCCTTAGGGTCACTCT
TGGAGAACCACAATGGATTCTTCTGTCTCTATACTTTCTGGGGTCTTTCCCTCCCTCTCATGAGAGTA
AAGTCCCCAGCAGATTTGAGTGAGTCAGGGGAAACACCCTCTTCTACTTACACTTGTCTCGCCACCCC
CCATTCTCTCCTTGGTTTCTTGGCCCTTACTGTGGCTGTACTTCAGATCATCAAAGCAAGAGAAAAGAAAT
GTGCTGAGTGTCTTCTCCCTGTGCCTAAGGGCCCTCATCCCAGCAGCCATATAGCCAGCAAATGGAG
GACTCCTAACTTCAGCACAAGTGGGTTTGGAGTGGCAGAGGCCATACTTTCTTGGCTCAGCTTCTCTCC
TTCTGCTGCCATCCTCTGCCTTAGAAGAAATGGGAACAGCTGTAGGGAAGGGACAAGCTGCAGAGG
AACAGTGATTTAAAGGCCCTTGGAGCTAGTCCCCTCAGCCTCACTTTGATCTGGCCTTAGCCTCAAAG
CTGCTGTACCTGCACATGGTCTGAGAACATCTCTTAGCCCTGTACCACACCTCTCTGCCTGTAGAACC
TGTGAGATCATTGCCTGTGGCCTCTGTACCCTGGAGCAAAGCCAGACCTTGGAAACCCAACTACTTTAG
TCCTAGAAAAGTTTCTGCTCAGAAGTGGAGTGGAAAGAGCCGTAAATGACGTTGGCCTTGGCCTGTTC
TTCTCCTCAGCGCCACACTCACATATTTCTGCTTTGGTTCTACAATCTGTTGTTGAGTCCCCCTCTC
TGCATCAGTCTTCCCTGGAAATGGAGGACCTGGTGCAGTACCTTCTCTGAGGCTTGGATTCCCCTTCT
TGGGTGTAATGTTGCTCCCTGATCCCATATCCCTGTGAAAATGACTTAGGAAAAACGTTTCTGAGCAAG
CATTCTGGCCCATGCACTGGCTGAAGTGGGCGACTGGAGGCGCCTCTGCAAGAAATGGAGTCAGGCGA
AGTCCACCTTGTACGCTGCTGTGCGCTTAGGAATAAAGTGAATCTCTTCCAACTGCAAAAAAAAAAAAAA AAAAAA

Restriction Sites: AscI-NotI



[View online >](#)

ACCN:	NM_019817
Insert Size:	534 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC058524 , AAH58524
RefSeq Size:	1686 bp
RefSeq ORF:	534 bp
Locus ID:	56447
UniProt ID:	P61924
Cytogenetics:	15 58.74 cM
Gene Summary:	The coatomer is a cytosolic protein complex that binds to dilysine motifs and reversibly associates with Golgi non-clathrin-coated vesicles, which further mediate biosynthetic protein transport from the ER, via the Golgi up to the trans Golgi network. Coatomer complex is required for budding from Golgi membranes, and is essential for the retrograde Golgi-to-ER transport of dilysine-tagged proteins. In mammals, the coatomer can only be recruited by membranes associated to ADP-ribosylation factors (ARFs), which are small GTP-binding proteins; the complex also influences the Golgi structural integrity, as well as the processing, activity, and endocytic recycling of LDL receptors (By similarity).[UniProtKB/Swiss-Prot Function]