

## Product datasheet for MC203537

### A730017C20Rik (NM\_173759) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	A730017C20Rik (NM_173759) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	A730017C20Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC075669

```
CAACACTGGCAAACCTTGTCTACTCACAGAGCCCGTGGCTTCTCAGAGGAACCTCTGTGAGCCGGAGAG
AAAAGGCGGCGAATCTGCATTATGCCTGCATACCAAGAAGACAGAGAGGAAGACACATGTACTATTGGGG
AGAGAGAAAAACCATTGCTACTCAGCGAATCAAGGGATTTGGTGTGGAGAACCCAGTCCTTCTGAGAG
TCTCCAGCATCCATGTCTTCAGTTATGAAGAACAACCCGTTGTATGGTGACATAAATTTGGAAGAAGCT
ATGGAAGAAAAGAAAAAGAGTCCTTCATGGACCATTGAGGAGTATGACAAGCATTCTGTGCATACAAATC
TCTCTGGACATCTGAAGGAAAACCCCAACGACCTGCGCTTTTGGTTGGGGGATACGTACACTCCTGGCTT
CGACACTCTGCTGAAAAAGAAGAAGAAACGCAACAAGCGATCAAATTTGTGCACATGGGTCTGATTTTA
CTCCTTGTGGCCTCCATTCTGGTTACCATAGTGACTCTCAGCACCATATTCTCTTAAGGAAGCTGCAGAC
CCATCTTTAAAGCTCTAAATGCCTTCTAAGTACTTGAAGAGAATCTGTCTGTGCATTGTAACGTCGG
CATTGCTAACTGTGTGACCAACACCATTGTCCTTTAGAATGTGAGCGTCGCTTGTGTTATGTTGGCCTC
TCTGCTTTGGTCTGTGAAACCTGTGTGTCATTCCAGTCTGGGTTTGTCTTTGTGAAGCAGGTATTCCAG
GGGCTCACAATATCAAAAGTCACTCCATACAGTTAACAATCAAGCTTTAGATTGGAAGCGGACTGCGCCA
AGTGCTGGCAACTGCTTGTCTTCATGCTTACCTGTTTTTACATCATTCTCCTTCTCCTCCTCCTCCTC
TCCTTCTCCTCCTCCTTCTCCTCCTTCTTTCTCCTCCTGCGCACTCGGGAAGCTGTAGCTCATCT
GCCCAACAGCTAGGGTGACCAGAAGGCTTCTGCTGAGCTACCCCTCAGTTTAGATCTGTTTTACCCTCCT
TTATCTCCGCTTAAGTTAAAGGCGGATAAAGTCGAAGGGAGACATTGTTCAAACCTATTGTAATTTCC
AGTGATCCTTTCTGAAACCTAAATTGTATGTTGTTGTTTTAGCAGAAAACAACTTGGTCAGCCACAGCT
GCATAATCAATTAAGTTACATCCAATCGTGCTTGGGGATTGATAGGAAGTGGCAAGTGTGACCTCTAGTG
GCAGTGGGTACACGTTTACCAGAAATAGGTGGAGGGGCACTGGAAGAGTCCAGGAATGTGCCACCC
TGAGCATTCAATTAAGTGAAGTACAGTCAAGCCGTAACAAGATTGTGAGAGAAGGGAAGTGGCTTAAGGATT
CTGTACAACATGGAGGGGAGAAGAGCGGTTGGCTCCTAACCTGTAGCTTTTTCTTCAGTGGTGATGTCGT
TTTTTTGTGATTGACAACCTTCTGGTGGGTGCAGATACCTGGGAACAAGCCCCAGGATACAAAGAACC
AAAGTGTCTGTCTTCTCCTGGACGTTCTTACCTAAGGCTCTCAGCGCCACTCCCAATGCTATTTAGTGT
GTGAAATAAGGCTTTTCTTCTTCTTCCCAAGAGTGCCGCTCTCGTGCCCTTGATAACGCCTTA
CTATCTCGCTGTATTGCTTATAACGATACTTGAAGTAAAGTATATAAAAAATGAGGGTTTACTGAAAAATA
AATAATTGTGTACAAAGCAAAAAAAAAAAAAAAAAAAAAA
```

Restriction Sites: Ascl-NotI



[View online >](#)

ACCN:	NM_173759
Insert Size:	456 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
RefSeq:	<a href="#">BC075669</a> , <a href="#">AAH75669</a>
RefSeq Size:	1787 bp
RefSeq ORF:	456 bp
Locus ID:	225583
UniProt ID:	<a href="#">Q8C4X7</a>
Cytogenetics:	18 D3